EIASON[®] hPTH

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EIASON® hPTH

Immunoenzymetric assay for the in vitro quantitative measurement of intact human Parathyroid Hormone (PTH) in serum and plasma


Instruction for Use

For in-vitro use only




Product of

IASON GmbH
Feldkirchner Straße 4
A – 8054 Graz-Seiersberg
Tel.: +43 (0)316 28 43 00
Fax: +43 (0)316 28 43 00-113
e-mail: order@iason.eu
www.iason.eu

Used IFU-Symbols

Symbol	English	Symbol	English
	In vitro diagnostic device		Conjugate
	Order number		Microtiterplate
	Product of		Wash solution
	Storage		Calibrators
	European Conformity		Control sera
	Expiry date		Substrate
	Batch code		Stop solution
	Microplate		Incubation Buffer

Intended use

For in-vitro use only.

The EIASON® hPTH kit is an immunoenzymetric assay for the in vitro quantitative measurement of intact human Parathyroid Hormone (PTH) in serum and plasma.

Summary

Biological Activity

Human parathyroid hormone (hPTH) is a major physiological regulator of phosphocalcic metabolism. hPTH increases serum calcium concentrations by its actions on kidney (enhancing tubular Ca^{++} reabsorption and phosphate excretion) and bone (stimulating osteoclastic activity and bone resorption). It indirectly affects intestinal absorption of Ca^{++} by stimulating renal 1α -hydroxylation of 25 hydroxyvitamin D. The release of PTH is controlled in a negative feedback loop by the serum concentration of Ca^{++} . PTH is synthesized in the chief cells of the parathyroid glands and secreted as an 84 amino acid molecule called "intact PTH", which is the main bioactive product. This molecule is degraded by proteolytic cleavage between amino acids 33-37 at peripheral sites to form biologically active aminoterminal fragments and biologically inactive carboxyl-terminal fragments. The carboxylterminal fragments are cleared only by glomerular filtration, while the bioactive intact PTH and amino-terminal fragments are also metabolically degraded in the liver and other tissues. The half-life of the carboxylterminal fragments increases dramatically in patients with renal failure. Thus, the measurement of intact PTH correlates best with the hormone production and biological activity.

Clinical application

The measurement of intact hPTH is used to establish the diagnosis of primary hyperparathyroidism by demonstrating elevated serum levels of bioactive PTH. It allows documenting the occurrence of secondary hyperparathyroidism in patients with Vit.D deficiency, intestinal malabsorption, or renal failure. In this last situation, the absence of interference with the inactive carboxyl-terminal fragments is especially valuable. The specificity and high sensitivity of the assay also allows differentiating clearly low serum PTH levels in hypoparathyroidism or in tumor-induced hypercalcaemia.

Assay principle

The EIASON® hPTH kit is a solid phase Enzyme Amplified Sensitivity Immunoassay performed on microtiterplates. Calibrators and samples react with the capture polyclonal antibodies (PAb, goat anti 1-34 PTH fragment) coated on microtiter well. After incubation, the excess of antigen is removed by washing. Then monoclonal antibodies (MAb, mouse anti 44-68 PTH fragment) labelled with horseradish peroxidase (HRP) are added. After an incubation period allowing the formation of a sandwich: coated PAb – human PTH – MAb – HRP, the microtiterplate is washed to remove unbound enzyme labelled antibody. Bound enzyme-labelled antibody is measured through a chromogenic reaction. The chromogenic solution (TMB) is added and incubated. The reaction is stopped with the addition of Stop Solution and the microtiterplate is then read at the appropriate wavelength. The amount of substrate turnover is determined colourimetrically by measuring the absorbance, which is proportional to the PTH concentration.

A calibration curve is plotted and PTH concentration in samples is determined by interpolation from the calibration curve. The use of the EASIA reader (linearity up to 3 OD units) and a sophisticated data reduction method (polychromatic data reduction) result in a high sensitivity in the low range and in an extended calibration range.

Warnings and precautions

The EIASON® hPTH kit is for in vitro diagnostic use only and is not for internal use in humans or animals. This product must be used strictly in accordance with the instructions set out in the Package Insert. IASON will not be held responsible for any loss or damage (except as required by statute) caused, arising out of non-compliance with the instructions provided.

CAUTION: this kit contains material of human and/or animal origin. Handle kit reagents as if capable of transmitting an infectious agent.

Source material from Human origin which is used in this kit was tested and found negative for HbsAG and HIV as well as for HCV antibodies. However, since there is no diagnostic procedure that excludes these agents with 100 percent certainty all components should be handled as potentially hazardous material.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. Disposal of kit reagents should be in accordance with local regulations.

Damaged test kit

In case of serious damage to the test kit or components, the company IASON must be notified in writing at least one week after receiving the kit. Severely damaged single components should not be used for the test run. They have to be stored until a final solution has been found.

Shelf life and storage of reagents

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the vial label, if kept at 2 to 8°C.
- Unused strips must be stored, at 2-8°C, in a sealed bag containing a desiccant until expiration date.
- After reconstitution, **CAL** and **CO** should be frozen immediately after use and kept at -20°C for three months. Avoid subsequent freeze-thaw cycles.
- The concentrated **WASH** is stable at room temperature until expiration date.
- Freshly prepared working wash solution should be used on the same day.
- After its first use, the **CONJ** is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

Specimen collection and preparation

- Blood samples should be promptly separated from the blood cells.
- Serum and plasma must be kept at 2 - 8°C.
- If the test is not run within 8 hours, storage in aliquots at -20°C is recommended. Avoid subsequent freeze thaw cycles.
- Prior to use, all samples should be at room temperature. It is recommended to vortex the samples before use.
- It is advisable to assay serum samples.
- Do not use haemolysed samples.

Materials provided

Allow all reagents to reach room temperature before use.

1. **MPL**: **Microtiterplate** with 96 anti PTH (polyclonal antibodies) coated breakable wells; ready to use.
2. **CAL0**: **Zero Calibrator** in human plasma and thymol; 1 vial lyophilized; add 3.0 mL aqua dest.
3. **CAL1**-**CAL5**: **Calibrators** in human plasma and thymol; 5 vials lyophilized; add 1.0 mL aqua dest.; exact concentration see bottles.
4. **CO1**; **CO2**: **2 controls** in human plasma with thymol; 2 vials lyophilized; add 1.0 mL aqua dest.
5. **WASH**: **Wash solution** NaCl-Tween20; 1 vial 25 mL; dilute 28 x with aqua dest. Prepare an adequate volume of working solution by adding 27 volumes of distilled water to 1 volume of **WASH** (28x). Use a magnetic stirrer to homogenize. Discard unused working wash solution at the end of the day.
6. **CONJ**: **Conjugate**; HRP labelled anti-PTH (monoclonal antibodies) in TRIS-maleate buffer with bovine serum albumin; thymol and sheep serum; 1 vial 11 mL; ready for use.
7. **INCBU**: **Incubation buffer** with EDTA, benzamidin and azide < 0.1%; 1 vial 6 mL; ready to use.
8. **SUB**: **Chromogenic TMB Solution** (Tetramethylbenzidine); 1 x 25 mL; ready to use;
9. **STOP**: **Stop solution**; HCl 1 N; 1 vial; 25 mL; ready to use.

Note:

1. Use **CAL0** for sample dilutions.
2. 1 pg of the calibrator preparation is equivalent to 1 pg of a synthetic PTH (1-84) from the Japanese Peptide Institute.

Materials required but not provided in the kit

- High quality aqua dest.
- Pipettes for delivery of: 50 µl, 100 µl, 200 µl, 1 mL, 2 mL and 3 mL (the use of accurate pipettes with disposable plastic tips is recommended)
- Vortex mixer
- Magnetic stirrer
- Horizontal microtiterplate shaker capable of 700 rpm ± 100 rpm
- Washer for microtiterplates
- Microtiterplate reader capable of reading at 450 nm, 490 nm and 650 nm (in case of polychromatic reading) or capable of reading at 450 nm and 650 nm (monochromatic reading)

Assay procedure

General remarks

- Do not use the kit or components beyond expiry date.
- Do not mix materials from different kit lots.
- Bring all the reagents to room temperature prior to use.
- Thoroughly mix all reagents and samples by gentle agitation or swirling.
- Perform **CAL**, **CO** and samples in duplicate. Vertical alignment is recommended.
- Use a clean plastic container to prepare the **WASH**.
- In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.
- For the dispensing of the **SUB** and the **STOP** avoid pipettes with metal parts.
- High precision pipettes or automated pipetting equipment will improve the precision.
- Respect the incubation times.
- To avoid drift, the time between pipetting of the first **CAL** and the last sample must be no longer than 60 minutes.
- Prepare a calibration curve for each run, do not use data from previous runs.
- Dispense the **SUB** within 15 minutes following the washing of the **MPL**.
- During incubation with **SUB**, avoid direct sunlight on the **MPL**.

Test procedure

1. Select the required number of strips for the run. The unused strips should be resealed in the bag with a desiccant and stored at 2-8°C.
2. Secure the strips into the holding frame.
3. Pipette 50 µL of **INCBU** into all wells.
4. Pipette 200 µL of each **CAL**, **CO** and sample into the appropriate wells.
5. Incubate for 2 hours at room temperature on a horizontal shaker set at 700 rpm ± 100 rpm.
6. Aspirate the liquid from each well.
7. Wash the plate 4 times by dispensing 400 µL of diluted **WASH** into each well and aspirating the content of each well.
8. Pipette 100 µL of **CONJ** into all the wells.
9. Incubate for 1 hour at room temperature on a horizontal shaker set at 700 rpm ± 100 rpm.
10. Aspirate the liquid from each well.
11. Wash the plate 4 times by dispensing 400 µL of diluted **WASH** into each well and aspirating the content of each well.
12. Pipette 100 µL of the **SUB** into each well within 15 minutes following the washing step.
13. Incubate the **MPL** for 30 minutes at room temperature on a horizontal shaker set at 700 rpm ± 100 rpm, avoid direct sunlight.
14. Pipette 200 µL of **STOP** into each well.
15. Read the absorbance at 450 nm and 490 nm (reference filter 630 nm or 650 nm) within 1 hour and calculate the results as next.

or fully automated on:

- **IASON® PersonalLab**
- **IASON® Quardette**
- **IASON® Gladiator**

Calculation of results

Polychromatic reading

In this case, the software will do the data processing.

1. The plate is first read at 450 nm against a reference filter set at 650 nm (or 630 nm).
2. A second reading is performed at 490 nm against the same reference filter.
3. The Software will drive the reader automatically and will integrate both readings into a polychromatic model. This technique can generate OD's up to 10.
4. The principle of polychromatic data processing is as follows:
Xi = OD at 450 nm
Yi = OD at 490 nm
Using a standard unweighted linear regression, the parameters A & B are calculated: $Y = A \cdot X - B$
If $X_i < 3$ OD units, then X calculated = X_i
If $X_i > 3$ OD units, then X calculated = $(Y_i - B) / A$
A 4 parameter logistic curve fitting is used to build up the calibration curve.
The PTH concentration in samples is determined by interpolation on the calibration curve.

Bichromatic reading

1. Read the plate at 450 nm against a reference filter set at 650 nm (or 630 nm).
2. Calculate the mean of duplicate determinations.
3. On semi-logarithmic or linear graph paper plot the OD values (ordinate) for each calibrator against the corresponding concentration of PTH (abscissa) and draw a calibration curve through the calibrator points by connecting the plotted points with straight lines.
4. Read the concentration for each control and sample by interpolation on the calibration curve.
5. Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a 4 parameter logistic function curve fitting is recommended.

Typical data

The following data are for illustration only and should never be used instead of the real time calibration curve.

EIASON® hPTH		OD _{units}
CAL	0.0 pg/mL	0.050
	22 pg/mL	0.149
	70 pg/mL	0.344
	224 pg/mL	0.999
	666 pg/mL	2.721
	1400 pg/mL	4.483

Expected Values

	pg/mL
Normal subjects	16 – 46
Hyperparathyroidism	106 - > 1000
Hypoparathyroidism	0 – 6.4

It is recommended that each laboratory determine a reference range for its own patient population. The range is based on 5% to 95% percentiles.

Quality control

The regular use of control samples at several analyte levels is advised to ensure day-to-day validity of results. Controls should be tested as unknowns. Quality Control charts should be maintained to follow the assay performance.

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid. In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or IASON GmbH directly.

Test characteristics

Detection limit

Twenty **CAL0** were assayed along with a set of other **CAL**. The detection limit, defined as the apparent concentration two standard deviations above the average OD at zero binding, was 2 pg/mL.

Specificity

Cross-reactive hormones or fragments were added to the zero calibrator, a high value calibrator (900 pg/mL) and a low value calibrator (100 pg/mL). The apparent PTH response was measured.

Cross-reactant	No significant interference up to
PTH 1-34 synthetic fragment, human	1000 pg/mL
PTH 44-68 synthetic fragment, human	20000 pg/mL
PTH 53-84 synthetic fragment, human	20000 pg/mL
PTH 73-84 synthetic fragment, human	100000 pg/mL
PTH-related protein 1-34 synthetic fragment, human	100000 pg/mL

Precision

Intra-Assay Precision				Inter-Assay Precision			
Serum	N	<X> ± SD (pg/mL)	CV (%)	Serum	N	<X> ± SD (pg/mL)	CV (%)
A	10	41.0 ± 0.5	1.1	A	20	45.7 ± 3.3	7.1
B	10	594 ± 12	2.0	B	20	381 ± 11.1	2.9

Accuracy

Recovery test

Sample	Added PTH (pg/mL)	Recovery PTH (pg/mL)	Recovery (%)
Serum	371	333	90
Heparin Plasma	371	347	93
EDTA Plasma	371	350	94

Dilution test

Sample	Dilution	Theoretical concent. (pg/mL)	Measured concent. (pg/mL)
Serum	1/1	-	955
	1/2	477	506
	1/4	239	229
	1/8	119	124
	1/16	60	63
	1/32	34	34

Samples were diluted with **CAL0**.

Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 60 minutes after the calibrators have been added to the coated wells.

	0´	15´	30´	45´	60´
S1	158	146	166	144	137
S2	77	72	70	67	72
S3	328	320	323	342	357
S4	260	250	250	258	251

Legal aspects***Reliability of results***

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact IASON.

Therapeutic consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under Reliability of Results. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point Therapeutic Consequences are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

Useful publications

1. HABENER J.F., and POTTE J.T., Jr. (1978) "Biosynthesis of parathyroid hormone". New Engl. J. Med., 299, 11:580 and 299, 12:635.
2. MARTIN K.J., HRUSKA K.A., FREITAG J.J., KLAH S. and SLOTOPOLSKY E. (1979) "The peripheral metabolism of parathyroid hormone". New Engl. J. Med., 301, 20:1092.
3. GOLTZMAN D., HENDERSON B. and LOVERIDGE N. "Cytochemical bioassay of PTH. (1980) Characteristics of the assay and analysis of circulating hormone forms". J. Clin. Invest., 65:1309.
4. POTTS J.T. Jr., KRONENBERG H.M., ROSENBLATT M. (1982) "Parathyroid hormone : Chemistry, biosynthesis and mode of action". Adv. Protein Chem., 323.
5. HACKENG W.H.L., LIPS P., NETELENBOS J.C. and LIPS C.J.M. (1986) "Clinical implication of estimation of intact parathyroid hormone (PTH) versus total immunoreactive PTH in normal subjects and hyperparathyroid patients". J. Clin. Endocrinol. Metab., 63:447.
6. BOUILLON R., COOPMANS W., DE GROOTE D.E.H., RADOUX D., ELIARD P.H. (1990) "Immunoradiometric assay of Parathyrin with polyclonal and monoclonal region specific antibodies". Clin. Chem., 36/2:271-276.

Pipetting scheme

Allow all reagents to reach room temperature before use

1. Pipetting	INCBU	50µL	
2. Pipetting	CAL0 – CAL5 200µL	CO1, CO2 200µL	sample 200µL
3. Incubation	2 hours at room temperature (18 – 25°C) on a shaker (700 rpm)		
4. Washing	Aspirate or decant and wash with 400µL WASH repeat 3 times, remove remaining liquid by hitting the plate against paper towel		
5. Pipetting	CONJ	100µL	
6. Incubation	1 hour at room temperature (18-25°C) on a shaker (700 rpm)		
7. Washing	Aspirate or decant and wash with 400µL WASH repeat 3 times, remove remaining liquid by hitting the plate against paper towel		
8. Pipetting	SUB	100µL	
9. Incubation	30 min at room temperature (18-25°C) on a shaker (700 rpm)		
10. Pipetting	STOP	200µL	
11. Reading	450nm (RF 630/650nm) and/or 490 nm (RF 630/650 nm) Calculation: 4-Parameter or Point To Point		

Expected Values

	pg/mL
Normal subjects	16 – 46
Hyperparathyroidism	106 - > 1000
Hypoparathyroidism	0 – 6.4