EIASON[®] anti IA-2

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EIASON® anti IA-2

Enzyme immunoassay for the quantitative determination of autoantibodies against
IA-2 in human serum


Instruction for use

For in-vitro-diagnostic use only



Product of

IASON GmbH
Feldkirchner Straße 4
A – 8054 Graz-Seiersberg
Tel.: +43 (0)316 28 43 00
Fax: +43 (0)316 28 43 00-113
Email: office@iason.eu
www.iason.eu

Used IFU-Symbols

Symbol	English	Symbol	English
	In vitro diagnostic device		Calibrators
	Order number		Control serum
	Product of		Streptavidin Peroxidase
	Storage		Dilution buffer for 
	Expiry date		Biotin
	European conformity		Buffer for reconstitution of 
	Batch code		Substrate
	Microplate		Stop solution
	Reaction enhancer		Concentrated wash solution

Intended use

For in-vitro-diagnostic use only.

The EIASON® anti IA-2 kit is an enzymeimmunoassay intended for the quantitative determination of autoantibodies against protein tyrosine phosphatase IA-2 in human serum. Autoantibodies to pancreatic beta cell antigens are important serological markers of type 1 diabetes mellitus (type 1 DM). The antigens recognised by these antibodies include insulin, glutamic acid decarboxylase (GAD₆₅ kDa isoform) and the islet cell antigen IA-2 or ICA-512.

Assay principle

In EIASON® anti IA-2 kit, IA-2 autoantibodies in patients' sera, calibrators and controls are allowed to interact with IA-2 coated onto ELISA plate wells. After a 16 - 20 hour incubation, the samples are discarded leaving IA-2 autoantibodies bound to the IA-2 coated wells. IA-2 Biotin is added in a 2nd incubation step where, through the ability of IA-2 autoantibodies to act divalently, a bridge is formed between the IA-2 immobilised on the plate and IA-2 Biotin. The amount of IA-2 Biotin is then determined in a third incubation step by the addition of Streptavidin Peroxidase, which binds specifically to Biotin. Excess, unbound Streptavidin Peroxidase is then washed away and addition of 3,3',5,5'-tetramethylbenzidine (TMB) results in formation of a blue colour. This reaction is stopped by addition of stop solution causing the well contents to turn from blue to yellow. The absorbance of the yellow reaction mixture at 405nm and 450nm is then read using an ELISA plate reader. A higher absorbance indicates the presence of IA-2 autoantibody in the test sample. Reading at 405nm allows quantitation of high absorbances (and should be used for concentrations of 200 U/mL or more). Low values (less than 30 U/mL) should be read off the 450nm calibrator curve. The measuring range is 15 – 4000 U/mL.

Warnings and precautions

The EIASON® anti IA-2 kit is for in vitro diagnostic use only and is not for internal use in humans or animals. This product must be used strictly in accordance with the instructions set out in the Package Insert. IASON will not be held responsible for any loss or damage (except as required by statute) caused, arising out of non-compliance with the instructions provided.

CAUTION: this kit contains material of human and/or animal origin. Handle kit reagents as if capable of transmitting an infectious agent.

Source material from Human origin which is used in this kit was tested and found negative for HbsAG and HIV as well as for HCV antibodies. However, since there is no diagnostic procedure that excludes these agents with 100 percent certainty all components should be handled as potentially hazardous material.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. Disposal of kit reagents should be in accordance with local regulations.

Safety considerations

SAPOD

Signal word: Warning

Hazard statement(s):

H317: May cause an allergic skin reaction

Precautionary statement(s):

P280: Wear protective gloves/protective clothing/ eye protection/face protection.

P302 + P352: IF ON SKIN: Wash with plenty of soap and water.

P333 + P313: If skin irritation or rash occurs: Get medical advice/attention.

P362 + P364: Take off contaminated clothing and wash it before reuse

SUB

Signal word: Danger

**Hazard statement:**

H360: May damage fertility or the unborn child

Precautionary statement(s):

P280: Wear protective gloves/protective clothing/eye protection/face protection

Damaged test kit

In case of serious damage to the test kit or components, the company IASON must be notified in writing at least one week after receiving the kit. Severely damaged single components should not be used for the test run. They have to be stored until a final solution has been found.

Shelf life and storage of reagents

This kit is stable until the stated expiry date if stored as specified. Upon receipt, store all reagents at 2-8°C.

Storage and preparation of serum samples

Sera to be analysed should be assayed soon after separation or stored, preferably in aliquots, at or below -20°C. 100 µL is sufficient for one assay (duplicate 50 µL determinations). Repeated freeze thawing or increases in storage temperature must be avoided. Do not use lipaemic or haemolysed serum samples. Do not use plasma in the assay. When required, thaw test sera at room temperature and mix gently to ensure homogeneity. Centrifuge serum prior to assay (preferably for 5 min at 10-15,000 rpm in a microfuge) to remove particulate matter. Please do not omit this centrifugation step if sera are cloudy or contain particulates.

Materials provided

1. **MPL** with wells coated with IA-2 (96 wells in total, 8 wells per strip).
Ensure stripwells are firmly fitted into frame provided. After opening return any unused wells to the original foil packet with desiccant provided and seal with adhesive tape. Place foil bag in the self-seal plastic bag and store at 2-8°C for up to 16 weeks.
2. **CAL 1** - **CAL 6** (0.7 mL each); ready to use: 0, 15, 60, 200, 400, 4000 U/mL corresponding to WHO reference preparation NIBSC 97/550.

3. **CO 1** (0.7 mL); ready to use. Refer to package insert for concentration.
4. **WASH** (125 mL); concentrated: Dilute 10 X with pure water before use. Store at 2–8°C up to expiry date.
5. **TURBO** reaction enhancer (4 mL, coloured red); ready to use.
6. **BIOBU** buffer for reconstituting **BIO** (2 x 15 mL, coloured blue); ready to use.
7. **BIO** (3 vials; lyophilised; reconstitute each vial with 5 mL of **BIOBU**; 1 vial for 4 strips). Use on day of reconstitution.
8. **PODBU** Diluent for diluting **SAPOD** (15 mL); ready to use.
9. **SAPOD** (Streptavidin-peroxidase; 0.7 mL); to be diluted 1:20 with **PODBU** (e.g. 0.5 mL + 9.5 mL). Store at 2-8°C for up to 20 weeks after reconstitution.
10. **SUB** Peroxidase substrate (tetramethyl benzidine; TMB; 15 mL ready to use).
11. **STOP** Stop solution (0.5 M sulphuric acid; 12 mL); ready to use.

Materials required but not provided in the kit

- Pipettes capable of dispensing 25 µL, 50 µL and 100 µL.
- Means of measuring out various volumes to reconstitute or dilute reagents.
- Pure water.
- ELISA Plate reader suitable for 96 well formats and capable of measuring at 450 nm and 405 nm.
- ELISA Plate shaker, capable of 500 shakes/min (not an orbital shaker).
- ELISA Plate cover.

Assay procedure

General remarks

1. All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
2. Once the test has been started, all steps should be completed without interruption.
3. Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination
4. Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
5. As a general rule the enzymatic reaction is linearly proportional to time and temperature.
6. Each assay run must include a standard curve and controls.

Test procedure

Allow all reagents reach room temperature (20-25°C) before use Do not reconstitute **BIO** until step 5.

All **CAL**, **CO** and samples, and should be run in duplicate and should be run concurrently so that all conditions of testing are the same.

1. Pipette 50 µL of **CAL 1** - **CAL 6**, **CO1** and test sera into the wells (in duplicate). Leave one well empty for blank.
2. Pipette 25 µL of **TURBO** into each well (except blank).
3. Cover the frame and shake the wells containing the various samples for 5 seconds on an ELISA plate shaker with 500 shakes per min.
4. Incubate the plate for 16–20 hours at 2–8 °C without shaking.
5. After the incubation, discard the samples by briskly inverting the frame of strip wells over a suitable receptacle. Wash 3 times with 300 µL of diluted **WASH** per well and each time tap the inverted wells gently on a clean dry absorbent surface to remove any droplets of wash buffer (not necessary when an automatic plate washer is used).
6. Pipette carefully 100 µL of reconstituted **BIO** into each well(except blank).
7. Incubate for 1 hour at **18-22°C** on an ELISA plate shaker (500 shakes per min). Incubation temperature **must not** exceed 22°C)
8. After the 1 hour incubation with **BIO**, wash 3 times with diluted **WASH** as described under point 5.
9. Pipette carefully 100 µL of diluted **SAPOD** into each well(except blank).
10. Incubate for 20 minutes at room temperature (20-25°C) on an ELISA plate shaker (500 shakes per min).
11. Wash 3 times with diluted **WASH** as described under point 5. If manual washing is being carried out use one additional wash step with pure water (to remove any foam) before finally tapping the inverted wells dry.
12. Pipette carefully 100 µL of **SUB** into each well (including blank).
13. Incubate in the dark for 20 minutes at 20-25°C without shaking.
14. Stop the substrate reaction by careful addition of 100 µL of **STOP** to each well (this will cause the blue colour to turn yellow) and shake the plate for about 5 seconds on a plate shaker to ensure uniformity of the solution in each well. It is most important to ensure that the substrate incubation time (i.e. time from addition of substrate to addition of stop solution) is the same for each well.
15. Measure the absorbance of each well at 450 nm and 405 nm for cover range filter (reference 620 – 650 nm) within 5 minutes after adding the **STOP**. Subtract the blank value.

or fully automated on:

1. **IASON® PersonalLab**
2. **IASON® Gladiator**

Calculation of results

- Draw a standard curve by plotting the absorbance of each standard against its concentration. Read off the values of the test samples.

- Alternative data reduction techniques may be employed but users should confirm that the selected curve fit is appropriate and gives acceptable results. 4PL (4 parameter logistics) or point-to-point fits are recommended.
- Samples with high IA-2 Ab concentrations can be 10x diluted in **CAL 1**. Please note that some sera will not dilute in a linear way.

Sample assay data


Typical results obtained with EIASON® anti IA-2 calibrators:

Calibrator U/mL	Absorbance	
	450 nm	405 nm
0	0.094	0.029
15	0.290	0.085
60	0.820	0.240
200	3.00	0.890
400	5.52	1.62
4000	8.33	2.45

This data is for illustration only and must not be used for the calculation of any sample result.

Absorbance reading at 405 nm can be converted to 450nm absorbance values by multiplying by the appropriate factor (3.0 in the case of IASON instrument).

Typical Calibration Curve

This sample calibration curve is for illustration only.

Expected Values

NIBSC 97/550	
	U/mL
Negative	<15
Positive	≥15

Each laboratory is recommended to determine ranges for their local population. The results alone should not be the only reason for any therapeutic consequences. The results should be correlated to other clinical observations and diagnostic tests.

Quality control

The regular use of control samples at several analyte levels is advised to ensure day-to-day validity of results. Controls should be tested as unknowns. Quality Control charts should be maintained to follow the assay performance.

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid. In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or IASON GmbH directly.

Test characteristics

Clinical specificity and sensitivity

In the DASP 2005 study the IASON IA-2 Ab ELISA kit showed 99% (n=100) specificity and 66% (n=50) sensitivity.

Lower detection limit

The kit negative control was assayed 20 times, and the mean and standard deviation calculated. The lower detection limit at +2 standard deviations was 0.3 U/mL

Precision

Inter-assay-precision (n=20)			Intra-assay-precision (n=25)		
Sample	Concentration [U/mL]	CV [%]	Sample	Concentration [U/mL]	CV [%]
A	143	3.5	A	142	3.6
B	33	6.6	B	36	51

Clinical accuracy

Analysis of sera from patients with autoimmune diseases other than type 1 DM indicated no interference from autoantibodies to thyroglobulin; thyroid peroxidase; acetylcholine receptor; TSH receptor or from rheumatoid factor.

Limitation of use

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice.

Any improper handling of samples or modification of this test might influence the results

Interference:

No interference was observed when samples were spiked with the following materials:

- Bilirubin up to 20 mg/dL
- Haemoglobin up to 5 mg/mL

Interference was observed with Intralipid at 1000 and at 3000 mg/dL.

Legal aspects**Reliability of results**

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact IASON.

Therapeutic consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under Reliability of Results. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point Therapeutic Consequences are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

Useful publications

1. S. Chen et al. Sensitive non isotopic assays for autoantibodies to IA-2 and to a combination of both IA-2 and GAD₆₅. Clinica Chimica Acta 2005 357:74-83
2. E. Nielsson et al. Calcium addition to EDTA plasma eliminates falsely positive results in the GADAb ELISA. Clinica chimica Acta 388 (2008) 130-134
3. K. Rahmati et al. A Comparison of Serum and EDTA Plasma in the Measurement of Glutamic Acid Decarboxylase Autoantibodies (GADA) and Autoantibodies to Islet Antigen-2 (IA-2A). Using the Radioimmunoassay (RIA) and Enzyme Linked Immunosorbent Assay (ELISA) Kits. Clin. Lab.2008 54:227-335
4. C. Törn et al. Diabetes Antibody Standardization Program: evaluation of assays for autoantibodies to glutamic acid decarboxylase and islet antigen-2. Diabetologia 2008 51:846-852
5. European patent 1448 993 B1, Chinese patent ZL02822274.1, Indian patent 226484 and related patents pending in other countries apply.

Patents

The following patents apply:

European patent EP 1 448 993 B1, Chinese patent ZL 02822274.1, Indian patent 226484, Japanese patent 5711449 and US patents US 8,129,132 B2, US 9,435,797 B2 and US 10,488,410 B2.

Pipetting scheme

Allow all reagents to reach room temperature (20-25°C) before use. Do not reconstitute **BIO** until step 6 below.

1. Pipetting	CAL 1 – CAL 6 50 µL	CO 1 50 µL	Samples 50 µL
2. Pipetting	TURBO	(except blank)	25 µL
3. Shaking	5 sec. on a shaker (500 shakes per min)		
4. Incubation	16-20 hours at 2-8°C without shaking		
5. Washing	wash 3x : aspirate or decant add 300µL of diluted WASH aspirate or decant and dry on an absorbent material*		
6. Pipetting	BIO	(except blank)	100µL
7. Incubation	1 hour at 18-22°C (500 shakes per min)		
8. Washing	wash 3x : see step 5*		
9. Pipetting	SAPOD	(except blank)	100µL
10. Incubation	20 min.at room temperature (20-25°C); on a shaker (500 shakes per min)		
11. Washing	wash 3x: see step 5		
12. Pipetting	SUB	(including blank)	100µL
13. Incubation	20 minutes at room temperature (20-25°C) in the dark without shaking		
14. Pipetting	STOP	(including blank)	100µL
15. Shaking	5 sec. on a shaker (500 shakes per min)		
16. Reading	450nm (RF 620nm) Overrange Filter: 405nm, Factor: 3 (dependent on photometer), reading within 5 min. Calculation: 4-parameter or point to point		

For details, see point assay procedure

Expected Values

NIBSC 97/550	
	U/mL
Negative	<15
Positive	≥15