EIASON[®] C-Peptide

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EIASON® C-Peptide

IVD

Enzyme immunoassay for the quantitative determination of C-Peptide
in human serum, plasma and urine

Kit instruction


For in-vitro use only

Product of


IASON GmbH

Feldkirchner Straße 4
A – 8054 Graz-Seiersberg
Tel.: +43 (0)316 28 43 00
Fax: +43 (0)316 28 43 00-4
Email: office@iason.eu
www.iason.eu

Used IFU-Symbols

Symbol	English	Symbol	English
	In vitro diagnostic device		Substrate
	Order number		Calibrators
	Product of		Stopping solution
	Storage		Wash buffer
	European Conformity		Batch code
	Expiry date		Sample diluent
	Enzym conjugate		Antiserum
	Enzyme complex		Microplate

Intended Use

For in-vitro use only.

The EIASON® C-Peptide ELISA is an enzyme immunoassay for the quantitative in vitro diagnostic measurement of C-Peptide in serum, plasma and urine.

Summary

Insulin is synthesized in the pancreatic beta cells as a 6000 MW component of an 86 amino acid polypeptide called proinsulin (1, 2, 3). Proinsulin is subsequently cleaved enzymatically, releasing insulin into the circulation along with a residual 3000 MW fragment called connection ("C") peptide, so-named because it connects A and B chains of insulin within the proinsulin molecule (1, 2, 3, 4). Human C-Peptide, a 31 amino acid residue peptide, has a molecular mass of approximately 3000 daltons.

C-Peptide has no metabolic function. However, since C-Peptide and insulin are secreted in equimolar amounts, the immunoassay of C-Peptide permits the quantitation of insulin secretion (4, 5, 6). This is the reason for the clinical interest of serum and urinary determinations of C-Peptide. Moreover, C-Peptide measurement has several advantages over immunoassays of insulin. The half-life of C-Peptide in the circulation is between two and five times longer than that of insulin (7). Therefore, C-Peptide levels are a more stable indicator of insulin secretion than the more rapidly changing levels of insulin. A very clear practical advantage of C-Peptide measurement arising from its relative metabolic inertness as compared to insulin is that C-Peptide levels in peripheral venous blood are about 5-6 times greater than insulin levels (3). Also, relative to an insulin assay, the C-Peptide assay's advantage is its ability to distinguish endogenous from injected insulin. Thus, low C-Peptide levels are to be expected when insulin is diminished (as in insulin-dependent diabetes) or suppressed (as a normal response to exogenous insulin), whereas elevated C-Peptide levels may result from the increased β -cell activity observed in insulinomas (3, 6, 9). C-Peptide has also been measured as an additional means for evaluating glucose tolerance and glibenclamide glucose tests (2, 3, 9, 10). C-Peptide levels are in many ways a better measurement of endogenous insulin secretion than peripheral insulin levels. C-Peptide may be measured in either blood or urine (9). With improved sensitive C-Peptide immunoassays, it is now possible to measure C-Peptide values at extremely low levels. The clinical indications for C-Peptide measurement include diagnosis of insulinoma and differentiation from factitious hypoglycemia, follow-up of pancreatectomy, and evaluation of viability of islet cell transplants (11, 12, 13). Recently, these indications have been dramatically expanded to permit evaluation of insulin dependence in maturity onset diabetes mellitus.

Clinical Indications for the EIASON® C-Peptide ELISA

- ◇ Assessment of residual β -cell function in diabetics under insulin therapy
- ◇ Detection and monitoring of the remission phase of type I diabetes
- ◇ Adjunct in the differential diagnosis between type I (insulin dependent) and type II (non-insulindependent) diabetes
- ◇ Diagnosis of insulin-induced factitious hypoglycemia.
- ◇ Contribution to the diagnosis of insulinoma (insulin suppression test)
- ◇ Prognostic index of fetal outcome in pregnant diabetic women
- ◇ Evaluation of insulin secretion in liver disease
- ◇ Monitoring of pancreatectomy

Assay principle

The EIASON® C-Peptide ELISA Kit is a solid phase enzyme-linked immunosorbent assay (ELISA), based on the principle of competitive binding. The microtiter wells are coated with anti-mouse antibodies, which bind a monoclonal antibody directed towards a unique antigenic site on the C-Peptide molecule. Endogenous C-Peptide of a patient sample competes with a C-Peptide-horseradish peroxidase conjugate for binding to the coated antibody. After incubation the unbound conjugate is washed off. The amount of bound peroxidase conjugate is inversely proportional to the concentration of C-Peptide in the sample. After addition of the substrate solution, the intensity of colour developed is inversely proportional to the concentration of C-Peptide in the patient sample.

Warnings and precautions

The EIASON® C-Peptide ELISA Kit is for in vitro diagnostic use only and is not for internal use in humans or animals. This product must be used strictly in accordance with the instructions set out in the Package Insert. IASON will not be held responsible for any loss or damage (except as required by statute) caused, arising out of non-compliance with the instructions provided.

CAUTION: this kit contains material of human and/or animal origin. Handle kit reagents as if capable of transmitting an infectious agent.

Source material from Human origin which is used in this kit was tested and found negative for HbsAG and HIV as well as for HCV antibodies. However, since there is no diagnostic procedure that excludes these agents with 100 percent certainty all components should be handled as potentially hazardous material.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. Disposal of kit reagents should be in accordance with local regulations.

Shelf life and storage of reagents

This kit is stable until the stated expiry date if stored as specified. Upon receipt, store all reagents at 2-8°C. Do not use reagents beyond this date.

Opened reagents must be stored at 2-8°C. Microtiter wells must be stored at 2-8°C. Once the foil bag has been opened, care should be taken to close it tightly again.

Reconstitute the lyophilized contents of the calibrator vials with 0.75 mL Aqua dest.

Note: The reconstituted calibrators are stable for 3 days at 2 °C to 8 °C. For longer storage freeze at -20°C.

The diluted Wash Solution is stable for 2 weeks at room temperature.

Specimen collection and preparation

Serum, plasma (EDTA-, heparin- or citrate plasma) or urine can be used in this assay.

Do not use haemolytic, icteric or lipaemic specimens.

Please note: Samples containing sodium azide should not be used in the assay.

Specimen collection

Serum:

Collect blood by venipuncture, allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

Plasma:

Whole blood should be collected into centrifuge tubes containing anti coagulant and centrifuged immediately after collection.

Urine:

The total volume of urine excreted during a 24 hour period should be collected and mixed in a single container.

Note: Specimens should be stored at 2-8°C during collection period and total volume collected should be recorded.

Specimen storage

Storage: Serum/Plasma

Specimens should be capped and may be stored for up to 24 hours at 2-8°C prior to assaying. Specimens held for a longer time should be frozen only once at -20°C prior to assay. When required, thaw test sera at room temperature and mix gently to ensure homogeneity.

Storage: Urine

Aliquot a well-mixed sample to be used in the assay. Centrifuge sample to clear. Urine samples may be stored for up to 36 hours at 2-8°C prior to assaying. Specimens held for a longer time should be frozen only once at -20°C prior to assay.

Specimen dilution

If in an initial assay, a specimen is found to contain more than the highest standard, the specimens can be diluted with **DIL** and reassayed as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

a) Dilution 1:10: 10 µL Serum + 90 µL **DIL** (mix thoroughly)

b) Dilution 1:100: 10 µL dilution a) 1:10 + 90 µL **DIL** (mix thoroughly).

Urine Dilution:

Prior to use dilute urine samples 1:20 with **DIL**

Urine samples

Prior to use dilute urine samples **1:20** with **DIL**

If the Sample Diluent included in the kit is insufficient, you can order additional **DIL** (40 mL vial) at IASON.

Materials provided

Allow all reagents 1-9 to reach room temperature before use.

10. **MPL** ELISA strip wells coated with anti mouse antibody. (96 wells in total, 8 wells per strip). Before opening the packet of strip wells, allow it to stand at room temperature (20-25°C) for at least 30 minutes. After opening, keep any unused wells in the original foil packet (reseal with adhesive tape) and in the self-seal plastic bag with the desiccant provided.
11. **CAL 0-5**; lyophilized: Reconstitute each vial with 0.75 mL Aqua. Dest. Stability: see "Storage of reagents". Concentrations: 0-16 ng/mL (see exact value on the vial label or on the QC-Datasheet). The calibrators are calibrated against IRR C-Peptide 84/510.
12. **DIL** Sample Diluent 1 vial, 3mL, ready for use
13. **AS** Antiserum, 1 vial, 7mL, ready for use
monoclonal mouse anti C-Peptide antibody

14. **CONJ** Enzyme Conjugate, 1 vial, 14 mL, ready to use biotinylated C-Peptide
15. **COMPL** Enzyme Complex, 1 vial, 14mL, ready for use contains horseradish peroxidase
16. **WASH** (30 mL, 40 x concentrated); dilute to 1.2 L with distilled water before use.
The diluted Wash Solution is stable for 2 weeks at room temperature
17. **SUB** Peroxidase substrate (tetramethyl benzidine; TMB; 14 mL ready to use).
18. **STOP** Stop solution (0.5 M sulphuric acid; 14 mL ready to use).

Note: Additional **DIL** for sample dilution is available upon request.

Materials required but not provided in the kit

- Pipettes capable of dispensing 50 µL and 100 µL
- Distilled water
- ELISA plate reader suitable for 96 well formats and capable of measuring absorbances at 450 nm.

Assay procedure

All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming. Once the test has been started, all steps should be completed without interruption. As a general rule the enzymatic reaction is linearly proportional to time and temperature. Calculate the number of individual **MPL** wells needed for the assay. Allow all the reagents supplied including the appropriate number of strips to reach room temperature, fit the number of strip wells required firmly into the frame provided. Controls and standard curve should always be included in each assay run.

1. Pipette 100 µL of **CAL 0-5**, controls and test sera into the wells (in duplicate).
2. Pipette carefully 50 µL of **ANTISER** into each well.
3. Pipette carefully 100 µL of **CONJ** into each well. Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step
4. Cover the frame and shake the wells containing the various samples and incubate for 60 minutes at 20-25° on an ELISA plate shaker.
5. After the incubation with **CONJ**, discard the **CONJ** by briskly inverting the wells over a suitable receptacle, wash 3 times with 400µL diluted **WASH** per well and each time tap the inverted wells gently on a clean dry absorbent surface to remove any droplets of wash buffer.
6. Pipette carefully 100 µL of **COMPL** into each well.
7. Cover the frame and shake the wells containing the various samples and incubate for 30 minutes at 20-25° on an ELISA plate shaker.

8. After the incubation with **COMPL**, discard the **COMPL** by briskly inverting the wells over a suitable receptacle, wash 3 times with 300µL diluted **WASH** per well and each time tap the inverted wells gently on a clean dry absorbent surface to remove any droplets of wash buffer.
9. Pipette carefully 100 µL of **SUB** into each well.
10. Incubate for 20 minutes at 20-25°C in the dark during which time a blue colour will develop.
11. Stop the substrate reaction by careful addition of 100 µL of **STOP** to each well (this will cause the blue colour to turn yellow) and shake the plate for about 5 seconds on a plate shaker to ensure uniformity of the solution in each well. It is most important to ensure that the substrate incubation time (i.e. time from addition of **SUB** to addition of **STOP**) is the same for each well.
12. Measure the absorbance of each well at 450 nm as reading filter (and 405 nm as overrange filter if possible) (reference 620 – 650 nm) within 10 minutes after adding the stop solution.

or fully automated on:

- **IASON® Quardette**
- **IASON® PersonalLab**
- **IASON® Gladiator**

Quality control

The regular use of control samples at several analyte levels is advised to ensure day-to-day validity of results. Controls should be tested as unknowns. Quality Control charts should be maintained to follow the assay performance.

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results. It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid. In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods. After checking the above mentioned items without finding any error contact your distributor or IASON directly.

Calculation of results

- Draw a standard curve by plotting the absorbance of each standard against its concentration. Read off the values of the test samples. Samples with concentrations higher than that of the highest standard have to be further diluted. For the calculation of the concentrations this dilution factor has to be taken into account.
- Alternative data reduction techniques may be employed but users should confirm that the selected curve fit is appropriate and gives acceptable results. 4PL (4 parameter logistics) or point-to-point fits are recommended.


Sample assay data

Example of typical results obtained with EIASON® C-Peptide calibrators:

Calibrator	Absorbance
ng / mL	450 nm
0	1.82
0.2	1.64
0.7	1.46
2	1.02
6	0.47
16	0.21

This data is for illustration only and must not be used for the calculation of any sample result.

Typical calibration curve

This sample calibration curve is for illustration only.

Expected values

	n	Mean ± 2SD
Serum (Post 12-hour fasting)	60	0.5 – 3.2 ng/mL
Urine		1 – 200 µg/day

Each laboratory is recommended to determine ranges for their local population. The results alone should not be the only reason for any therapeutic consequences. The results should be correlated to other clinical observations and diagnostic tests.

Performance characteristics

Assay dynamic range

The range of the assay is between 0.06 – 16 ng/mL.

Specificity to antibodies (cross-reactivity)

The cross-reactivity of intact or split-Proinsulin is clinically not significant.

Analytical sensitivity: 0.064 ng/mL

The analytical sensitivity was calculated from the mean plus two standard deviations of twenty (20) replicate analyses of **CAL 0**

Reproducibility

Intra-assay-precision (n = 20)			Inter-assay-precision (n = 12)		
Sample	Concentration (ng/mL)	CV [%]	Sample	Concentration (ng/mL)	CV [%]
1	0.48	6.54	1	0.42	9.33
2	2.30	6.70	2	2.05	9.92
3	3.86	5.13	3	4.23	8.38

Recovery

Samples have been spiked by adding C-Peptide solutions with known concentrations in a 1:1 ratio. The % Recovery has been calculated by multiplication of the ratio of the measurements and the expected values with 100.

Serum sample	Endogenous C-Peptide ng/mL	Added concentration 1:1 (v/v) (ng/mL)	Measured conc. (ng/mL)	Expected conc. (ng/mL)	Recovery (%)
1	5.36	0.0	5.36		
		16.0	10.31	10.68	96.6
		6.0	5.57	5.68	98.0
		2.0	3.63	3.68	98.7
		0.7	3.08	3.03	101.8
2	9.7	16.0	9.7	12.85	97.2
		6.0	12.49	7.85	104.8
		2.0	8.23	5.85	87.9
		0.7	5.15	5.20	87.2
3	12.12	16.0	12.12	14.06	110.4
		6.0	15.52	9.06	107.3
		2.0	9.72	7.06	103.4
		0.7	7.30	6.41	88.1
		0.7	5.65		

Urine sample	Endogenous C-Peptide ng/mL	Added concentration 1:1 (v/v) (ng/mL)	Measured conc. (ng/mL)	Expected conc. (ng/mL)	Recovery (%)
1	4.1	16.0	2.1	10.1	107.9
		6.0	10.	5.1	109.2
		2.0	5.57	3.05	85,2
2	2.1	16.0	1.01	9.01	102.1
		6.0	9.2	4.0	100.5
		2.0	4.03	2.01	109.5
3	5.0	16.0	2.5	10.5	96.2
		6.0	10.1	5.5	96.4
		2.0	5.3	3.5	108.6
		2.0	3.8		

Linearity

Serum sample	Dilution	Measured conc. (ng/m)	Expected conc. (ng/m)	Recovery (%)
A	None	6.1	6.1	
	1:2	3.25	3.05	106.7
	1:4	1.61	1.52	105.3
	1:8	0.84	0.76	110.6
	1:16	0.41	0.38	107.6
B	None	9.90	9.90	
	1:2	5.59	4.95	112.8
	1:4	2.48	2.48	100.3
	1:8	1.29	1.24	104.0
	1:16	0.69	0.62	111.8
C	None	13.25	13.25	
	1:2	6.97	6.62	105.1
	1:4	3.22	3.31	97.1
	1:8	1.7	1.66	102.8
	1:16	0.85	0.83	103.1

Urine sample	Dilution	Measured conc. (ng/m)	Expected conc. (ng/m)	Recovery (%)
A	None	8.7	8.7	
	1:2	4.29	4.335	98.6
	1:4	2.01	2.18	92.4
	1:8	1.09	1.09	100.2
B	None	9.2	9.2	
	1:2	4.7	4.6	102.2
	1:4	2.25	2.3	97.8
	1:8	1.12	1.15	97.5
C	None	13.9	13.9	
	1:2	6.6	6.95	95.0
	1:4	3.3	3.48	95.0
	1:8	1.8	1.74	103.6

Limitation of use

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice. Any improper handling of samples or modification of this test might influence the results.

Interferences:

Haemoglobin (up to 4 mg/m), Bilirubin (up to 0.5 mg/m) and Triglyceride (up to 30 mg/m) have no influence on the assay results.

Drug Interferences:

Until today no substances (drugs) are known to us, which have an influence to the measurement of C-Peptide in a sample.

No hook effect as observed in this test.

Legal aspects

Reliability of results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact IASON GmbH.

Therapeutic consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point Reliability of Results. Any laboratory result is only a part of the total clinical picture of a patient. Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived. The test result itself should never be the sole determinant for deriving any therapeutic consequences.

Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point Therapeutical Consequences are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

Useful publications

1. Ashby, J. and Frier, B.: Circulating C-Peptide: Measurement and Clinical Applications. *Annals of Clinical Biochemistry*. 18:125, 1981
2. Beischer, W.: Proinsulin and C-Peptide in Humans. *Hormones in Normal and Abnormal Human Tissues*. Volume 3K, Fotherby and Pal, S., ed. (Berlin: Walter DeGruyter). pp. 1-43, 1983
3. Beyer, J., Krause V., Cordes V.: C-Peptide: Its Biogenesis, Structure, Determination and Clinical Significance. *Giornale Italiano di Chimica Clinica* 4 Supp. 9:22, 1979
4. Bonger, A. and Garcia-Webb, P.: C-Peptide Measurement: Methods and Clinical Utility. *CRC Critical Reviews in Clinical Laboratory Sciences*. 19:297, 1984.
5. Blix, P. Boddie-Wills, C., Landau, R., Rochman, H. Rubenstein, A.: Urinary C-Peptide: An Indicator of Beta-Cell Secretion under Different Metabolic Conditions. *Journal of Clinical Endocrinology and Metabolism*. 54:574, 1982.
6. Rendell, M.: C-Peptide Levels as a Criterion in Treatment of Maturity-Onset Diabetes. *Journal of Clinical Endocrinology and Metabolism*. 57 (6): 1198, 1983
7. Horwitz, D., et al.: Proinsulin, Insulin and C-Peptide concentrations in Human Portal and Peripheral Blood. *Journal of Clinical Investigation*. 55:1278, 1975
8. Horwitz, D., Kurzuya, H., Rubenstein, A.: Circulating Serum C-Peptide. *The New England Journal of Medicine*. 295:207, 1976
9. Rendell, M.: The Expanding Clinical Use of C-Peptide, Radioimmunoassay. *Acta Diabetologica Latina*. 20:105, 1983
10. Heding, L. and Rasmussen, S.: Human C-Peptide in Normal and Diabetic Subjects. *Diabetologica*. 11:201, 1975
11. Canivet, B., Harter, M., Viot, G., Balgrac, N., Krebs, B.: Residual β -Cell Function in Insulin-Dependent Diabetes: Evaluation by Circadian Determination of C-Peptide Immuno reactivity. *Journal of Endocrinological Investigation*. 3:107, 1980.
12. Starr, J., Horwitz, D., Rubenstein, A., Mako, M.: Insulin, Proinsulin and C-Peptide. *Methods of Hormone Radioimmunoassay* 2nd Ed., Academic Press Inc., 1979
13. Rubenstein, A., Kuruya, H., Horwitz, D.: Clinical Significance of Circulating C-Peptide in Diabetes Mellitus and Hypoglycemic Disorders. *Archives of Internal Medicine*. Vol. 137:625, May 1977.
14. Yalow, R., Berson, S.: Introduction and General Considerations. *Principles of Competitive Protein Binding Assays*. Ch. 2, Eds. Odell, W. and Daugheday, W., J.B. Lippincott Co., Philadelphia, 1971

Pipetting scheme

Allow all reagents to reach room temperature (20-25°C)

1. Pipetting	CAL 0-5 100 µL	Controls 100 µL	Samples 100µL
2. Pipetting	AS		50µL
3. Pipetting	ECONJ		100µL
	Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step		
4. Incubation	1 hours at room temperature on a shaker (20-25°C)		
5. Washing	wash 3 x : aspirate or decant add 400 µL WASH aspirate or decant and dry on an absorbent material		
6. Pipetting	COMPL		100µL
7. Incubation	30 minutes at room temperature on a shaker (20-25°C)		
8. Washing	wash 3 x : aspirate or decant add 400 µL WASH aspirate or decant and dry on an absorbent material		
9. Pipetting	SUB		100µL
10. Incubation	20 min at room temperature (20-25°C)		
11. Pipetting	STOP		100µL
12. Reading	450nm (RF 620 – 650nm) Optional Overage Filter: 405nm, Factor: 3 (dependent on photometer), reading within 10 min. Calculation: 4-parameter or point to point		

Expected values

NIBSC 84/510	
Serum (Post 12-hour fasting)	0.5 – 3.2 ng/mL
Urine	1 – 200 µg/day