EIASON[®] AChRAb

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EIASON® AChRAb

Enzymeimmunoassay for the quantitative determination of acetylcholine receptor autoantibodies (AChRAb) in human serum

Kit instruction

For in-vitro use only

Product of

IASON GmbH

Feldkirchner Straße 4


A – 8054 Graz-Seiersberg

Tel.: +43 (0)316 28 43 00


Fax: +43 (0)316 28 43 00-4

Email: office@iason.eu

www.iason.eu

Used IFU-Symbols

Symbol	English	Symbol	English
	In vitro diagnostic device	CAL1 – CAL4	Calibrators
REF	Order number	CO1 CO2	Positive controls
	Product of	NC	Negative control
	Storage	BIO	MAB Biotin
	European Conformity	BIOBU	Buffer for BIO
	Expiry date	SAPOD	Streptavidin peroxidase
	Lot number	PODBU	Buffer for SAPOD
AChRFOET	Foetal type AChR	SUB	Diluent for SAPOD
AChRBU	Buffer for AChRFOET	STOP	Stop Solution
AChRADULT	Adult type AChR	WASH	Concentrated wash solution
MPL	Microplate		

Intended use

For in-vitro use only.

The EIASON® AChRAb kit is an enzymeimmunoassay intended for the quantitative determination of autoantibodies against acetylcholine receptor (AChR) in human serum and plasma, responsible for failure of the neuromuscular junction in myasthenia gravis.

Assay principle

EIASON® AChRAb depends on the ability of AChRAb in human serum to bind to similar sites on the AChR as various monoclonal antibodies such as MAb1 (coated on ELISA plate wells) and/or MAb2 and/or MAb3 (which are labelled with biotin). In the absence of AChRAb a complex is formed between MAb1 coated on the plate wells, the AChR and MAb2 and MAb3 biotin. MAb2 and MAb3 Biotin bound are then detected by addition of streptavidin peroxidase, substrate and stop solution. The higher the concentration of AChRAb in the test serum, the greater the inhibition of MAb-biotin binding.

Warnings and precautions

The EIASON® AChRAb kit is for in vitro diagnostic use only and is not for internal use in humans or animals. This product must be used strictly in accordance with the instructions set out in the Package Insert. IASON will not be held responsible for any loss or damage (except as required by statute) caused, arising out of non-compliance with the instructions provided.

CAUTION: this kit contains material of human and/or animal origin. Handle kit reagents as if capable of transmitting an infectious agent.

Source material from Human origin which is used in this kit was tested and found negative for HbsAG and HIV as well as for HCV antibodies. However, since there is no diagnostic procedure that excludes these agents with 100 percent certainty all components should be handled as potentially hazardous material.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. Disposal of kit reagents should be in accordance with local regulations.

Shelf life and storage of reagents

This kit is stable until the stated expiry date if stored as specified. Upon receipt, store all reagents at 2-8°C.

Storage and preparation of serum samples

Sera to be analysed should be assayed soon after separation or stored, preferably in aliquots, at or below -20°C. 100µL is sufficient for one assay (duplicate 50µL determinations). Repeated freezing and thawing or increases in storage temperature must be avoided. Incorrect storage of serum samples can lead to loss of autoantibody activity.

Do not use lipaemic or grossly haemolysed serum samples.

Studies in which EDTA, citrate and heparin plasma samples were spiked with AChRAb positive sera showed minor changes in signal compared with spiked serum from the same donor. In particular OD450 values with spiked EDTA, citrate and heparin plasmas were 83% - 122% of spiked serum (20 samples with serum concentrations ranging from 0.28 nmol/L – 18 nmol/L) or 69% - 165% in terms of nmol/L.

When required, bring test sera to room temperature and mix gently to ensure homogeneity. Centrifuge serum prior to assay (preferably for 5 min at 10-15000 rpm in a microfuge) to remove particulate matter. Please do not omit this centrifugation step if sera cloudy or contain particulates.

Materials provided

Allow all reagents 1-14 to reach room temperature before use.

1. **AChRFOET** Foetal Type AChR; 3 vials; lyophilised. Reconstitution with 0.7 mL **AChRBU**. When more than one vial is used, pool the vials and mix gently before use. Pool the vials when more than one vial is required and then use immediately to reconstitute **AChRADULT**.
2. **AChRADULT** Adult Type AChR; 3 x 0.5 mL each; lyophilised. Reconstitution with 0.5 mL **AChRFOET** solution to give a mixture of **AChRFOET** and **AChRADULT**. When more than one vial is used, pool the vials and mix gently before use. Use up to 6 hours after reconstitution if stored at 2-8°C.²
3. **AChRBU** Buffer for **AChRFOET**; 1 x 5 mL; ready for use.
4. **MPL** Microplate coated with MAb1; 12 breakapart strips of 8 wells (96 in total) in a frame and sealed in a foil bag. Allow foil bag to stand at room temperature (20-25°C) for 30 minutes before opening. Ensure wells are firmly fitted in the frame provided. After opening return any unused wells to the original foil bag and seal with adhesive tape. Then place foil bag in the selfseal plastic bag with desiccant provided, and store at 2-8°C for up to kit expiry date.
5. **CAL1** – **CAL4** Calibrators; 0.5; 1.0; 6.5; 20 nmol/L, toxin bound; 0.7 mL each, ready for use.
6. **CO1**, **CO2** Positive Controls; 0.7 mL each; ready for use; see label for concentration range.
7. **NC** Negative Control; 3.0 mL; ready for use.
8. **BIO** MAb Biotin (MAb2 and MAb3); 3 vials; lyophilised; reconstitute each vial with the volume of **BIOBU** shown on the vial label. Mix gently, and leave to stand at room temperature (20–25 °C) for 5 minutes before use. When more than one vial is used, pool the vials and mix gently before use. Store at 2 – 8°C for up to kit expiry date after dilution.
9. **BIOBU** Buffer for **BIO**; 15 mL; ready for use.

² The absorbance at 450 nm will be 10-15% lower when reconstituted receptors have been stored for 6 hours at 2-8C.

10. **SAPOD** Streptavidin Peroxidase, 0.7 mL; concentrated; dilute 1 in 20 with **PODBU**. For example, 0.5 mL **SAPOD**+ 9.5 mL **PODBU**. Store at 2 – 8°C for up to 16 weeks after dilution.
11. **PODBU** Diluent for **SAPOD**; 15 mL; ready for use.
12. **SUB** Peroxidase Substrate TMB; ready for use; 15 mL.
13. **WASH** Concentrated Wash Solution; 100 mL; concentrated; Dilute 1 in 10 with pure water before use. Use up to expiry of kit.
14. **STOP** Stop Solution; 10 mL; ready for use.

Materials required but not provided in the kit

- Pipettes capable of dispensing 25 µL, 50 µL and 100 µL
- Eppendorf tubes
- Distilled water
- ELISA plate reader suitable for 96 well formats and capable of measuring absorbances at 450 and 405 nm
- ELISA plate shaker (capable of 500 shakes/min), not an orbital shaker

Assay procedure

Calculate the number of individual **MPL** wells needed for the assay. Allow all the reagents supplied including the appropriate number of strips to reach room temperature, fit the number of strip wells required firmly into the frame provided. Kit controls should always be included in each assay run.

Please note that the optical density depends on incubation time and temperature. Therefore, it is necessary to bring all reagents to work-ready state before the start of the assay; the caps should be opened and all required wells should be in strip holder. Only such preparation will ensure equal elapsed time for each pipetting step without interruption.

DAY ONE

1. Pipette 100 µL **NC** **CAL1** - **CAL4**, **CO1**, **CO2** and test sera into 1.5 mL Eppendorf tubes.
2. Pipette 25 µL of **AChRFOET** and **AChRADULT** mixture into each Eppendorf tube and seal the tubes. Centrifuge for 10 s at 10 – 15,000 rpm. Vortex gently and incubate overnight (16 – 20 h) at 2 – 8°C.

DAY TWO

1. Gently mix each tube using a vortex mixer. Pipette 50 µL of each sample-AChR mixture into the ELISA plate wells. Cover the wells and incubate at room temperature (20 – 25°C) on an ELISA plate shaker (500 shakes per min) for one hour.
2. After the 1-hour incubation discard the samples by briskly inverting the frame of strip wells over a suitable receptacle. Wash 3 times with 300 µL

- WASH** per well and each time tap the inverted wells gently on a clean dry absorbent surface to remove any droplets of wash buffer.
3. Pipette carefully 50 µL of reconstituted **BIO** into each well.
 4. Cover the plate and incubate for 1 hour at 20-25°C on an ELISA plate shaker (500 shakes per min).
 5. After the 1 hour incubation with **BIO**, discard the **BIO** by briskly inverting the wells over a suitable receptacle, wash 3 times with **WASH** as described under point 4.
 6. Pipette carefully 100 µL of **SAPOD** into each well.
 7. Cover the plate and incubate for 30 minutes at room temperature (20-25°C) on an ELISA plate shaker (500 shakes per min).
 8. Discard the **SAPOD** by briskly inverting the wells over a suitable receptacle; wash 3 times with **WASH** as described under point 4. For manual washing, wash once more with pure water to remove any foam.
 9. to remove any foam. followed by once with pure water to remove any foam.
 10. Pipette carefully 100 µL of peroxidase **SUB** (TMB) into each well.
 11. Incubate for 30 minutes at room temperature (20-25°C) in the dark during which time a blue colour will develop.
 12. Stop the substrate reaction by careful addition of 50 µL of **STOP** to each well (this will cause the blue colour to turn yellow) and shake the plate for about 5 seconds on a plate shaker to ensure uniformity of the solution in each well. It is most important to ensure that the substrate incubation time (i.e. time from addition of **SUB** to addition of **STOP**) is the same for each well.
 13. Measure the absorbance of each well at 450 nm and 405 nm for overrange filter (reference 620 – 650 nm) within 30 minutes after adding the stop solution.

or fully automated on:

- **IASON® Gladiator**
- **IASON® Quardette**
- **IASON® Personal Lab**

Quality control

The regular use of control samples at several analyte levels is advised to ensure day-to-day validity of results. Controls should be tested as unknowns. Quality Control charts should be maintained to follow the assay performance.

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid. In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or IASON GmbH directly.

Calculation of results

- Draw a standard curve by plotting the absorbance of each standard against its concentration. Read off the values of the test samples.
- Alternative data reduction techniques may be employed but users should confirm that the selected curve fit is appropriate and gives acceptable results. 4PL (4 parameter logistics) or point-to-point fits are recommended.
- Samples with high AChRAb concentrations can be diluted in EIASON® AChRAb [NC]. For example 10 µL of sample plus 90 µL of [NC] to give a 10 x dilution. Some sera will not dilute in a linear way and we suggest that the dilution giving a value closest to 50 % inhibition is used for calculation of AChRAb concentration.


Sample assay data

Typical results obtained with EIASON® AChRAb calibrators:

Calibrators	450 nm	nmol/L
NC	1.970	0.2
CO1	0.469	7.5
CO2	1.124	1.6
CAL1	1.616	0.5
CAL2	1.329	1.0
CAL3	0.524	6.5
CAL4	0.144	20

This data is for illustration only and must not be used for the calculation of any sample result.

Typical calibration curve

This sample calibration curve is for illustration only. Results can also be expressed as inhibition (%) of AChR binding calculated using the formula;

$$100 \times \left(1 - \frac{\text{test sample absorbance at 450 nm}}{\text{negative control (G) absorbance at 450 nm}} \right)$$

This % inhibition value can then be converted to nmol/L toxin bound using the formula:

$$0.2 \times 2^{(0.067 \times \% \text{ inhibition of test sample})}$$

This formula has been established empirically using a comparison of AChRab measurements by the IASON ELISA and RIA methods. Close agreement between nmol/L values obtained in the AChRab ELISA using the calibration curve and using this formula should not be expected in the case of all individual sera.

Sample	Abs. 450nm	% Inhibition	Calculated nmol/L
NC	1.97	0	0.2
CO1	0.46	76.2	6.9
CO2	1.12	42.9	1.5

Expected Values

	nmol/L
Negative	< 0.45
Positive	≥ 0.45

Each laboratory is recommended to determine ranges for their local population.

Clinical specificity:

402 individual healthy blood donors were assayed in the EIASON® AChRAb ELISA. 401 (99.8 %) were identified as being negative for AChR autoantibodies. One sample was positive and gave a value of 20% inhibition (0.54 nmol/L from the calibration curve, 0.52 nmol/L calculated).

Clinical sensitivity:

Samples from 83 patients diagnosed with myasthenia gravis were assayed in the the EIASON® AChRAb ELISA. 76 (92%) were identified as being positive for AChR autoantibodies.

Lower detection limit

The negative control was assayed 20 times and the mean and standard deviation calculated. The lower detection limit at 2 standard deviations was 0.25 nmol/L

Precision

Inter assay precision (n = 20)					Intra assay precision (n = 24)				
Sample	% Inhibition	CV (%)	nmol/L	VK(%)	Probe	% Inhibition	CV (%)	nmol/L	VK(%)
1	76.4	3.3	7.7	8.7	1	90.8	0.6	13.5	2.5
2	52.4	6.7	2.0	11.1	2	45.9	2.4	1.7	5.2
3	27.3	9.4	0.62	9.4	3	25.9	7.1	0.67	8.5

Clinical accuracy

Analysis of 107 sera from patients with autoimmune diseases other than myasthenia gravis indicated no interference from autoantibodies to thyroglobulin (n=10), thyroid peroxidase (n=11), dsDNA (n=9), TSH receptor (n=40), glutamic acid decarboxylase (n=10), 21-hydroxylase (n=10), or from rheumatoid factor (n=27). Two other samples gave values of 28% (0.74 nmol/L) and 44% (1.5 nmol/L) inhibition and were from a patient with Systemic Lupus Erythematosus (dsDNA Ab positive) respectively. These samples were assayed in the IASON AChRAb RIA kit and were positive (values of 1.3 and 1.5 nmol/L respectively). In addition two samples from patients with rheumatoid arthritis (rheumatoid factor positive) were positive in the IASON AChRAb ELISA and gave values of 24% (0.77 nmol/L) and 19% (0.61 nmol/L) inhibition. The first of these samples was also positive in the IASON AChRAb RIA (5.3 nmol/L)

Assay specificity

There are no interferences concerning the following substances:

- Intralipide up to 3000 mg/dl
- Bilirubin up to 20 mg/dl
- Haemoglobin up to 500 mg/dl

LEGAL ASPECTS

Reliability of results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact IASON.

Therapeutic consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under Reliability of Results. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point Therapeutic Consequences are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

References

R. Hewer et al. A sensitive non-isotopic assay for acetylcholine receptor autoantibodies. Clinica Chimica Acta 2006 **364**: 159 - 166

This page is left blank intentionally

Pipetting scheme

Allow all reagents to reach room temperature before use

1. Pipetting	NC CAL1 – CAL4 , CO1 , CO2 , Samples in 1.5 mL Eppendorf Tubes 100 µL
2. Pipetting	AChRFOET and AChRADULT Mixture 25 µL Centrifuge for 10 s at 10-15.000 g Vortex gently Incubate overnight 16-20 hrs at 2-8°C
3. Pipetting	50 µL of each sample-AchR mixture on MPL
4. Incubation	1 h at room temperature (20-25°C) with shaking (500 rpm)
5. Washing	wash 3x : aspirate or decant add 300 µL WASH aspirate or decant and dry on an absorbent material
6. Pipetting	BIO 50 µL
7. Incubation	1 h at room temperature (20-25°C) with shaking (500 rpm)
8. Washing	wash 3x : aspirate or decant add 300 µL WASH aspirate or decant and dry on an absorbent material
9. Pipetting	SAPOD 100 µL
10. Incubation	30 min at room temperature (20-25°C)
11. Washing	wash 3x : aspirate or decant add 300 µL WASH aspirate or decant and dry on an absorbent material
12. Pipetting	SUB 100 µL
13. Incubation	30 min at room temperature (20-25°C)
14. Pipetting	STOP 50 µL
15. Reading	at 450 nm / 4Parameter or Point to Point within 30 min.

Expected Values

	nmol/L
Negative	< 0.45
Positive	≥ 0.45