# EIASON<sup>®</sup> Aquaporin-4 Ab V2

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# EIASON® Aquaporin-4 Ab V2

Enzyme immunoassay for the quantitative determination of autoantibodies against  
AQP4 in human serum

## Instruction for use

For in-vitro-diagnostic use only

Product of

**IASON** GmbH

Feldkirchner Straße 4

A – 8054 Graz-Seiersberg


Tel.: +43 (0)316 28 43 00

Fax: +43 (0)316 28 43 00-113


Email: [office@iason.eu](mailto:office@iason.eu)

[www.iason.eu](http://www.iason.eu)

REF E07-050-96-2

## Used IFU-Symbols

Symbol	English	Symbol	English
	In vitro diagnostic device		Positive Controls
	Order number		Negative Control
	Product of		Biotin
	Expiry date		Buffer for reconstitution of 
	Storage		Streptavidin Peroxidase
	European Conformity		Diluent for 
	Batch code		Concentrated Wash Solution
	Microplate		Substrate
	Calibrators		Stop solution

## Intended use

*For in-vitro-diagnostic use only.*

The EIASON® Aquaporin-4 Ab V2 kit is an enzymeimmunoassay intended for use by professional persons only, for the quantitative determination of autoantibodies against AQP4 in human serum. Neuromyelitis optica (NMO), also known as Devic's syndrome, is an immune-mediated neurologic disease that involves the spinal cord

and optic nerves. It can be considered to be a disorder distinct from multiple sclerosis (MS). A serum immunoglobulin G autoantibody (NMO-IgG) has been shown to be a specific marker for NMO and the water channel aquaporin 4 (AQP4) has been identified as the antigen for NMO IgG. Measurement of AQP4 autoantibodies can be of considerable value in distinguishing NMO from MS when full clinical features may not be apparent and early intervention may prevent or delay disability.

### **Assay principle**

In the EIASON® Aquaporin-4 Ab V2 Kit, AQP4 antibodies in patients' sera, calibrators and controls are allowed to interact with AQP4 coated onto ELISA plate wells and liquid phase biotinylated AQP4. After incubation at room temperature for 2 hours with shaking, the well contents are discarded. AQP4 Ab bound to the AQP4 coated on the well will also interact with AQP4-Biotin through the ability of AQP4 Ab in the samples to act divalently leaving AQP4-Biotin bound to the well via an AQP4 Ab bridge. The amount of AQP4-Biotin bound is then determined in a second incubation step involving addition of streptavidin peroxidase (SA-POD), which binds specifically to biotin. Excess, unbound streptavidin peroxidase is then washed away and addition of the peroxidase substrate, 3,3',5,5'- tetramethylbenzidine (TMB), results in formation of a blue colour. This reaction is stopped by the addition of a stop solution, causing the well contents to turn yellow. The absorbance of the yellow reaction mixture at 450nm and 405nm is then read using an ELISA plate reader. A higher absorbance indicates the presence of AQP4 autoantibody in the test sample. Reading at 405nm allows quantitation of high absorbances. It is recommended that values below 10 U/mL should be measured at 450nm. If it is possible to read at only one wavelength 405nm may be used. The measuring range is 3 – 80 U/mL (arbitrary IASON units).

### **Warnings and precautions**

The EIASON® Aquaporin-4 Ab V2 kit is for in vitro diagnostic use only and is not for internal use in humans or animals. This product must be used strictly in accordance with the instructions set out in the Package Insert. IASON GmbH will not be held responsible for any loss or damage (except as required by statute) caused, arising out of non-compliance with the instructions provided.

**CAUTION:** this kit contains material of human and/or animal origin. Handle kit reagents as if capable of transmitting an infectious agent.

Source material from Human origin which is used in this kit was tested and found negative for HbsAG and HIV as well as for HCV antibodies. However, since there is no diagnostic procedure that excludes these agents with 100 percent certainty all components should be handled as potentially hazardous material.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. Disposal of kit reagents should be in accordance with local regulations.

## Safety considerations

SAPOD

**Signal word:** Warning

**Hazard statement(s):**

H317: May cause an allergic skin reaction

**Precautionary statement(s):**

P280: Wear protective gloves/protective clothing/ eye protection/face protection.

P302 + P352: IF ON SKIN: Wash with plenty of soap and water.

P333 + P313: If skin irritation or rash occurs: Get medical advice/attention.

P362 + P364: Take off contaminated clothing and wash it before reuse

## Shelf life and storage of reagents

This kit is stable until the stated expiry date if stored as specified. Upon receipt, store all reagents at 2-8°C. Do not use reagents beyond this date.

## Storage and preparation of serum samples

Sera to be analyzed should be assayed soon after separation or stored, preferably in aliquots, at or below –20°C. 100 µL is sufficient for one assay (duplicate 50 µL determinations). Repeated freeze thawing or increases in storage temperature should be avoided. Do not use lipaemic or haemolysed samples. Studies in which EDTA, citrate and heparin plasma samples were spiked with AQP4 Ab positive sera showed minor changes in signal compared with spiked serum from the same donor. In particular OD<sub>450</sub> values with spiked EDTA, citrate and heparin plasmas were 79% - 128% of spiked serum (15 samples with serum concentrations ranging from 2.6 U/mL – 30 U/mL) or 87% - 130% in terms of U/mL.

When required, thaw test sera at room temperature and mix gently to ensure homogeneity. Centrifuge serum prior to assay (preferably for 5 min at 10 - 15,000 rpm in a microfuge) to remove particulate matter. Please do not omit this centrifugation step if sera are cloudy or contain particulates.

## Materials provided

1. **MPL** Microplate with wells coated with AQP4 (96 wells in total, 8 wells per strip). Before opening the packet of strip wells, allow it to stand at room temperature (20-25°C) for at least 30 minutes. After opening, keep any unused wells in the original foil packet (reseal with adhesive tape) and in the self-seal plastic bag with the desiccant provided. Store at 2-8°C and use within 4 months. A frame for holding the wells during assays is also provided.
2. **CAL 1** - **CAL 5** Calibrators; 0.7 mL each; ready to use; 1.5; 5; 20; 40; 80 U/mL, arbitrary IASON units.

3. **CO1** **CO2** Positive Controls; 0.7 mL each; ready to use. Refer to package insert for concentration.
4. **NC** Negative Control; 0.7 mL; ready to use.
5. **WASH** Concentrated wash solution; 120 mL; 10 x concentrated; dilute 1: 10 with distilled water before use; store at 2 – 8°C up to kit expiry date.
6. **BIOBU** AQP4 Biotin reconstitution buffer for reconstituting **BIO**; 1 x 10 mL; ready to use.
7. **BIO** AQP4 Biotin; 3 vials; lyophilised; reconstitute immediately before use each vial with 1.5 mL of **BIOBU**; when more than one vial is to be used, pool the vials and mix gently before use.
8. **PODBU** Diluent for diluting **SAPOD**; 15 mL; ready to use.
9. **SAPOD** Streptavidin Peroxidase; 0.8 mL; to be diluted 1: 20 with **PODBU** (e.g. 0.5 mL + 9.5 mL). Store at 2-8°C for up to 4 months after reconstitution.
10. **SUB** Peroxidase substrate; Tetramethyl benzidine; TMB; 15 mL ready to use.
11. **STOP** Stop solution; 14 ml; ready to use.

### Materials required but not provided in the kit

- Pipettes capable of dispensing 25 µL, 50 µL and 100 µL.
- Means of measuring out various volumes to reconstitute or dilute reagents.
- Pure water.
- ELISA Plate reader suitable for 96 well formats and capable of measuring at 450 nm and 405 nm.
- ELISA Plate shaker, capable of 500 shakes/min (not an orbital shaker).
- ELISA Plate cover.

### Assay procedure

#### General remarks

1. All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
2. Once the test has been started, all steps should be completed without interruption.
3. Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination
4. Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
5. As a general rule the enzymatic reaction is linearly proportional to time and temperature.
6. Each assay run must include a standard curve and controls.

### Test procedure

Allow all reagents to reach room temperature (20-25°C) before use.

All **CAL**, **CO** and samples, and should be run in duplicate and should be run concurrently so that all conditions of testing are the same.

1. Pipette 50 µL of **CAL 1** - **CAL5**, **NC**, **CO1** **CO2** and test sera into the wells (in duplicate). Leave one well empty for blank.
2. Pipette 25 µL of **BIO** into each well (except blank).
3. Cover the frame and shake the wells for 2 hours at room temperature (20 -25°C) on an ELISA plate shaker (500 shakes per min).
4. After the incubation, discard the samples by briskly inverting the frame of strip wells over a suitable receptacle. Wash 3 times with 300 µL diluted **WASH** per well and each time tap the inverted wells gently on a clean dry absorbent surface to remove any droplets of wash buffer.
5. Pipette carefully 100 µL of **SAPOD** into each well (except blank).
6. Cover the plate and incubate for 20 min at room temperature (20-25°C) on an ELISA plate shaker (500 shakes per min).
7. Repeat wash step 4.
8. Pipette carefully 100 µL of **SUB** into each well (including blank).
9. Incubate for 20 min in the dark at room temperature (20-25°C) without shaking.
10. Stop the substrate reaction by careful addition of 100 µl of **STOP** to each well (this will cause the blue colour to turn yellow) and shake the plate for about 5 seconds on a plate shaker to ensure uniformity of the solution in each well. It is most important to ensure that the substrate incubation time (i.e. time from addition of substrate to addition of stop solution) is the same for each well.
11. Measure the absorbance of each well at 450 nm and 405 nm for overrange filter (reference 620 – 650 nm) within 10 minutes after adding the **STOP**. Subtract the blank value.

or fully automated on:

- **IASON® PersonalLab**
- **IASON® Gladiator**

### Calculation of results

A calibration curve can be established by plotting calibrator concentration on the x-axis (log scale) against the absorbance of the calibrators on the y-axis (linear scale). The AQP4 autoantibody concentrations in patient sera can then be read off the calibration curve. Other data reduction systems can be used. The **NC** can be assigned a value of 0.15 U/mL to assist in computer processing of assay results. Samples with AQP4 autoantibody concentrations above 80 U/mL can be diluted (e.g.

10 x and/or 100 x) in AQP4 antibody negative serum. Some sera will not dilute in a linear way.

**Sample assay data**


Typical results obtained with EIASON® Aquaporin-4 Ab V2 calibrators:

Calibrators U/mL	Absorbance	
	450 nm	405 nm
1.5	0.085	0.024
5	0.354	0.108
20	1.577	0.468
40	3.802	1.118
80	6.588	1.938

This data is for illustration only and must not be used for the calculation of any sample result.

Absorbance readings at 450 nm can be converted to 405 nm absorbances by multiplying by the appropriate factor (3.0 in the case of IASON instruments).

**Typical calibration curve**

This sample calibration curve is for illustration only.

### Expected Values

	U/mL
Negative	<3
Positive	≥3

Each laboratory is recommended to determine ranges for their local population. The results alone should not be the only reason for any therapeutic consequences. The results should be correlated to other clinical observations and diagnostic tests.

### Quality control

The regular use of control samples at several analyte levels is advised to ensure day-to-day validity of results. Controls should be tested as unknowns. Quality Control charts should be maintained to follow the assay performance.

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid. In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or IASON GmbH directly.

### Test characteristics

#### *Clinical specificity*

Sera from 358 individual healthy blood donors were tested in the EIASON® Aquaporin-4 Ab V2 Kit. 356 (99%) sera were identified as being negative for AQP4 Ab

#### *Clinical sensitivity*

Of 62 sera from patients with NMO or NMO spectrum disorder (NMOSD) 48 (77%) were positive for AQP4 Ab.

**Lower detection limit**

The **NC** was assayed 20 times and the mean and standard deviation calculated. The lower detection limit at 2 standard deviations was 0.17 U/mL.

**Precision**

Inter-assay-precision (n=20)			Intra-assay-precision (n=25)		
Sample	Concentration [U/mL]	CV [%]	Sample	Concentration [U/mL]	CV [%]
1	5.0	15.6	1	3.9	7.7
2	13.3	10.5	2	7.0	8.6
3	35	7.9	3	28	3.2
4	59	7.5	4	58	3.1

**Clinical accuracy**

Analysis of 205 sera from patients with autoimmune diseases other than neuromyelitis optica spectrum disorders (NMOSD) indicated no interference from autoantibodies to the TSH receptor (n=110), glutamic acid decarboxylase (n=26), 21-hydroxylase (n=12), the acetylcholine receptor (n=10), thyroid peroxidase (n=15), thyroglobulin (n=10), IA-2 (n=7) or from rheumatoid factor (n=15) in the EIASON® Aquaporin-4 Ab V2 Kit.

**Limitation of use**

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice.

Any improper handling of samples or modification of this test might influence the results.

**Interference**

No interference was observed when samples were spiked with the following materials:

- Intralipid at 1000 and 3000 mg/dL
- Bilirubin up to 20 mg/dL

Interference was seen from haemoglobin at 500 mg/dL.

## **Legal aspects**

### ***Reliability of results***

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact IASON GmbH.

### ***Therapeutic consequences***

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under Reliability of Results. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

### ***Liability***

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point Therapeutic Consequences are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

## Literature

1. V.A. Lennon et al. A serum autoantibody marker of neuromyelitis optica: distinction from multiple sclerosis. *Lancet* 2004 364(9451): 2106-2112
2. V.A. Lennon et al. IgG Marker of optic-spinal multiple sclerosis binds to the aquaporin-4 water channel. *The journey of Experimental Medicine* 2005 202: 473 – 477
3. B.G. Weinshenker et al. Neuromyelitis optica IgG predicts relapse after longitudinally extensive transverse myelitis. *Annals of Neurology* 2006 59: 566 - 569
4. N. Isobe et al. Quantitative assays for anti-aquaporin-4 antibody with subclass analysis in neuromyelitis optica. *Multiple Sclerosis Journal* 2012 18: 1541 – 1551
5. S. Jarius et al. Testing for antibodies to human aquaporin-4 by ELISA: Sensitivity, specificity and direct comparison with immunohistochemistry. *Journal of the Neurological Sciences* 2012 320: 32 – 37

## Patents

The following patents apply:

European patent EP 1 700 120 B1, US patents US 7,101,679 B2, US 7,947,254 B2 and US 8,889,102 B2 , Chinese patent ZL200480040851.3 and Japanese patent 4538464.

## Pipetting scheme

Allow all reagents and sample to reach room temperature (20 – 25°C) before use.

1. Pipetting	CAL1 – CAL5 50µL	NC 50µL	CO1 – CO2 50µL	Samples 50µL
2. Pipetting	BIO	(except blank)		25µL
3. Incubation	Cover plate 2 hours at room temperature (20-25°C) with shaking (500 shakes per min)			
4. Washing	wash 3x : aspirate or decant add 300µL of diluted WASH aspirate or decant and dry on an absorbent material*			
5. Pipetting	SAPOD	(except blank)		100µL
6. Incubation	Cover plate 20 min at room temperature (20-25°C) with shaking (500 shakes per min)			
7. Washing	see point 4. *			
8. Pipetting	SUB	(including blank)		100µL
9. Incubation	20 min at room temperature (20-25°C) <b>in the dark</b> without shaking			
10. Pipetting	STOP	(including blank)		100µL
11. Reading	450nm (RF 620nm) Overrange Filter: 405nm, Factor: 3 (dependent on photometer), reading within 10 min. Calculation: 4-parameter or point to point			

For details, see point assay procedure

### Expected values

	U/mL
Negative	<3
Positive	≥3