# EIASON<sup>®</sup> FSME (TBE) IgG

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# EIASON<sup>®</sup> FSME (TBE) IgG

Enzyme immunoassay for the quantitative determination of IgG-class antibodies against TBE-Virus in human serum or plasma (citrate, heparin)


## Kit instruction

For in-vitro use only



Product of

**IASON GmbH**  
Feldkirchner Straße 4  
A – 8054 Graz-Seiersberg  
Tel.: +43 (0)316 28 43 00  
Fax: +43 (0)316 28 43 00-113  
Email: [order@iason.eu](mailto:order@iason.eu)  
[www.iason.eu](http://www.iason.eu)

## Used IFU-Symbols

Symbol	English	Symbol	English
	In vitro diagnostic device		Microplate
	Order number		Sample diluent
	Product of		Calibrators
	Storage		Substrate solution
	European Conformity		Enzym conjugate
	Expiry date		Stop solution
	Batch code		Washing buffer

## Intended use

*For in-vitro use only.*

The EIASON® FSME (TBE) IgG ELISA is intended for the quantitative determination of IgG class antibodies against TBE virus in human serum or plasma (citrate, heparin).

## Summary

Tick-borne encephalitis (TBE) virus is a flavivirus of the family *Togaviridae*. It is an enveloped single-stranded RNA virus with cubic icosahedral symmetry and ranges in size from 20-80nm in diameter.

Three subtypes can be distinguished which show only little differences in their structural proteins.

TBE virus is mainly transmitted by ticks. The degree of contamination of ticks (and thus humans) in central Europe increases from west to east, and anybody may be affected. Specific antibody development yields a life-long immunity.

TBE is the most important tick-transmitted disease of man -beside Lyme disease, which is caused by the spirochete *Borrelia burgdorferi*. The clinical course of the disease depends on the immune status of the infected persons. A high virus production in the primary infected tissues is required for the passage of the blood-brain barrier and the resulting severe manifestations in the central nervous system.

Species	Disease	Symptoms (e.g.)	Transmission route
TBE-Virus	Tick-borne encephalitis  CEE (Central European Encephalitis)	Phase 1: unspecific flu-like symptoms (mild fever, headache, muscle pain, joint pain, gastrointestinal complaints) Phase 2: high fever, development of meningitis and / or encephalitis	By tick bites ( <i>Ixodes ricinus</i> , western subtype; <i>Ixodes persulcatus</i> , eastern subtype). Rarely by infected (nonpasteurized) milk.

The presence of pathogen or infection may be identified by:

- Serology: e.g. ELISA, Neutralization, Hemagglutination Inhibition, Complement Fixation

### Assay principle

The quantitative immunoenzymatic determination of specific antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a second washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA microwell plate reader.

### Warnings and precautions

The EIASON® FSME (TBE) IgG ELISA Kit is for in vitro diagnostic use only and is not for internal use in humans or animals. This product must be used strictly in accordance with the instructions set out in the Package Insert. IASON will not be held responsible for any loss or damage (except as required by statute) caused, arising out of non-compliance with the instructions provided.

**CAUTION:** this kit contains material of human and/or animal origin. Handle kit reagents as if capable of transmitting an infectious agent.

Source material from Human origin which is used in this kit was tested and found negative for HbsAG and HIV as well as for HCV antibodies. However, since there is no diagnostic procedure that excludes these agents with 100 percent certainty all components should be handled as potentially hazardous material.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. Disposal of kit reagents should be in accordance with local regulations.

### Damaged test kit

In case of serious damage to the test kit or components, the company IASON must be notified in writing at least one week after receiving the kit. Severely damaged single components should not be used for the test run.

### Shelf life and storage of reagents

This kit is stable until the stated expiry date if stored as specified. Upon receipt, store all reagents at 2-8°C. Do not use reagents beyond this date.

### Storage and preparation of serum samples

Use human serum or plasma (citrate, heparin) samples with this assay. If the assay is performed within 5 days after sample collection, the specimen should be kept at 2 - 8°C; otherwise they should be aliquoted and stored deep-frozen (-70 – -20°C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing. Heat inactivation of samples is not recommended.

### Sample dilution

Before assaying, all samples should be diluted 1+100 with **DIL**. Dispense 10 µL sample and 1 mL **DIL** into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

### Materials provided

*Allow all reagents to reach room temperature before use.*

1. **MPL TBE-Virus Coated Wells (IgG)**: 12 breakapart 8-well snap-off strips coated with TBE-Virus antigen; in resealable aluminium foil. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2-8 °C.
2. **DIL IgG Sample Diluent**: 1 bottle containing 100 mL of phosphate buffer (10mM) for sample dilution; containing anti human-IgG; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap.
3. **STOP Stop Solution**: 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
4. **WASH Washing Solution (20 x conc.)**: 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2 for washing the wells; white cap. Dilute washing solution 1+19; e.g. 10 mL washing solution + 190 mL distilled water. The diluted buffer will keep for 5 days if stored at room temperature. Crystals in the solution disappear by warming up to 37 °C in a water bath. Mix well before dilution.
5. **CONJ TBE-Virus anti-IgG Conjugate**: 1 bottle containing 20 mL of peroxidase labelled antibody to human IgG in phosphate buffer; coloured blue; 10 mM; ready to use; black cap.
6. **SUB TMB Substrate Solution**: 1 bottle containing 15 mL 3,3',5,5'-tetra-methyl-benzidine (TMB), < 0.1 %; ready to use; yellow cap; < 5% NMP. The reagent is ready to use and has to be stored at 2- 8 °C, away from the light. The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.
7. **CAL0-CAL4 TBE Virus IgG Standards**: 5 vials, each containing 2 mL standard (human serum or plasma); coloured yellow; ready to use:  
**CAL 0**: 0 IAU/mL; blue cap

CAL 1:	50	IAU/mL; green cap	
CAL 2:	130	IAU/mL; yellow cap	[IAU = IASON arbitrary units]
CAL 3:	200	IAU/mL; red cap	
CAL 4:	300	IAU/mL; white cap	

For potential hazardous substances please check the safety data sheet

### Materials supplied

- 1 Cover foil
- 1 Instruction for use (IFU)
- 1 Plate layout

### Materials required but not provided in the kit

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450/620nm (or IASON fully automated instrument)
- Incubator 37°C
- Manual or automatic equipment for rinsing wells
- Pipettes to deliver volumes between 10 and 1000 µl
- Vortex tube mixer
- Deionised or (freshly) distilled water
- Disposable tubes
- Timer

### Assay procedure

#### General remarks

1. All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
2. Once the test has been started, all steps should be completed without interruption.
3. Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination
4. Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
5. As a general rule the enzymatic reaction is linearly proportional to time and temperature.
6. Each assay run must include a standard curve and controls.

#### Test procedure

Please read the instruction for use carefully before performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three to five and the volume of Washing Buffer from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 12. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to  $37 \pm 1$  °C.

1. Dispense 100 µL **CAL0**-**CAL4** and diluted samples into their respective wells. Leave well A1 for substrate blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour ± 5 min at 37±1°C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of **WASH**. Avoid overflows from the reaction wells. The soak time between each wash cycle should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!

*Note: Washing is important! Insufficient washing results in poor precision and false results.*

5. Dispense 100 µL **CONJ** into all wells except for the blank well (e.g. A1). Cover with foil.
6. **Incubate for 30 min at room temperature. Do not expose to direct sunlight.**
7. Repeat step 4.
8. Dispense 100 µL **SUB** into all wells.
9. **Incubate for exactly 15 min at room temperature in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL **STOP** into all wells in the same order and at the same rate as for the **SUB** thereby a colour change from blue to yellow occurs.
11. Measure the absorbance of the specimen at 450/620nm within 30 min after addition of the **STOP**.

or fully automated on:

- **IASON<sup>®</sup> Quardette**
- **IASON<sup>®</sup> PersonalLab**
- **IASON<sup>®</sup> Gladiator**

## Measurement

Adjust the ELISA Microwell Plate Reader **to zero** using the **substrate blank**.

If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!

**Measure the absorbance** of all wells at **450 nm** and record the absorbance values for **CAL1**-**CAL4** and patient sample in the plate layout.

Dual wavelength reading using 620 nm as reference wavelength is recommended.

Where applicable calculate the **mean absorbance values** of all duplicates.

## Run validation criteria

In order for an assay to be considered valid, the following criteria must be met:

<b>Substrate BLANK</b>	Absorbance < <b>0.100</b>
<b>CAL0</b>	Absorbance < <b>0.200</b>
<b>CAL1</b>	Absorbance > <b>0.050</b>
<b>CAL2</b>	Absorbance > <b>CAL1</b>

$\frac{\text{CAL3}}{\text{CAL4}}$  Absorbance >  $\frac{\text{CAL2}}{1.000}$   
 Absorbance > 1.000

$\text{CAL0} < \text{CAL1} < \text{CAL2} < \text{CAL3} < \text{CAL4}$

If these criteria are not met, the test is not valid and must be repeated.

### Calculation of results

In order to obtain quantitative results in IAU/mL plot the (mean) absorbance values of the 5 Standards  $\text{CAL0}$ - $\text{CAL4}$  on (linear/linear) graph paper in a system of coordinates against their corresponding concentrations (0, 50, 130, 200 and 300 IAU/mL) and draw a standard calibration curve (absorbance values on the vertical y-axis, concentrations on the horizontal x-axis).

Read results from this standard curve employing the (mean) absorbance values of each patient sample.

For the calculation of the standard-curve mathematical Point to Point function should be used.

### Results interpretation

	IAU/mL	
Grey zone	55 - 110	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks. If the result is unequivocal again the sample is judged as negative. <b>After vaccination:</b> This may be a case of seroconversion. Continue with basic immunisation or booster. Repeat test within 2-4 weeks. A non-specific reaction is not excluded
Negativ	< 55	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely. <b>After vaccination:</b> No seroconversion after vaccination. This may be the case after the first vaccination during basic immunisation. Patients which show low or no response.
Postiv	> 110	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine). <b>After vaccination:</b> This is a case of seroconversion. Check anamnestic data and if necessary complete basic immunization or give a booster.

Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data. In immunocompromised patients and newborns serological data only have restricted value.

**Antibody isotypes and state of infection**

---

<b>Serology</b>	<b>Significance</b>
IgM	Characteristic of the primary antibody response High IgM titer with low IgG titer: → suggests a current or very recent infection Rare: → persisting IgM
IgG	Characteristic of the secondary antibody response May persist for several years High IgG titer with low IgM titer: → may indicate a past infection

---

**Quality control**

The regular use of control samples at several analyte levels is advised to ensure day-to-day validity of results. Controls should be tested as unknowns. Quality Control charts should be maintained to follow the assay performance.

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid. In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or IASON GmbH directly.

**Test characteristics**

The results refer to the groups of samples investigated; these are not guaranteed specifications.

***Precision***

Intraassay Precision (n=24)			Interassay Precision (n=12)		
Probe	Mean [OD]	CV [%]	Probe	Mean [IAU]	CV [%]
1	1.411	3.51	1	26.55	9.18
2	1.991	2.23	2	15.38	12.99
3	1.913	8.93	3	253.81	10.48

***Diagnostic specificity***

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 100.0% (95% confidence interval: 93.84% - 100.0%).

***Diagnostic sensitivity***

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 100.0% (95% confidence interval: 96.31% - 100.0%).

***Analytical sensitivity***

The analytical sensitivity (according to CLSI EP17-A) is defined as the apparent concentration of the analyte that can be distinguished from the zero calibrator. It is 1.78 IAU/mL.

***Interferences***

Interferences with hemolytic, lipemic or icteric samples are not observed up to a concentration of 10 mg/mL hemoglobin, 5 mg/mL triglycerides and 0.5 mg/mL bilirubin.

***Cross reactivity***

Cross reactivity with other flaviviruses cannot be excluded and should be taken into account for result interpretation.

***Measurement range***

The measurement range is 1.78 IAU/mL – 300 IAU/mL.

***Limitation of use***

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

## **Legal aspects**

### ***Reliability of results***

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact IASON.

### ***Therapeutic consequences***

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under Reliability of Results. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

### ***Liability***

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point Therapeutic Consequences are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

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## Pipetting scheme

Allow all reagents to reach room temperature before use

Dilute all samples 1+100 with **DIL**.  
Optional: Leave A1 for substrate-Blank

1. Pipetting	<b>CAL 0-4</b> <b>100 µL</b>	Diluted samples <b>100 µL</b>
2. Inkubation	1 hour at 37°C	
3. Washing	wash 3 x: aspirate or decant add <b>300 µL</b> of diluted <b>WASH</b> aspirate or decant and dry on an absorbent material	
4. Pipetting	<b>CONJ</b> (except blank)	<b>100 µL</b>
5. Inkubation	30 min at room temperature (20-25°C)	
6. Washing	wash 3 x: See point 3.	
7. Pipetting	<b>SUB</b>	<b>100 µL</b>
8. Inkubation	15 min. at room temperature (20-25°C) in the dark	
9. Washing	<b>STOP</b>	<b>100 µL</b>
10. Messen	450nm (RF 620nm) Overrange Filter: 405nm, Factor 3 (depends on photometer), reading within 30 min. Calculation: 4 parameter or point to point	

\*Please increase the volume of diluted **WASH** to 350 µL if necessary

### Expected values

	IAU/mL
Negative	< 55
Grey zone	55 – 110
Positive	>110