EIASON[®] VZV IgG

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Literatur

1. Croen, Kenneth D.; Straus, Stephen E. (1991): Varicella-zoster virus latency. In *Annual review of microbiology* 45, pp. 265–282. DOI: 10.1146/annurev.mi.45.100191.001405.
2. Dlugosch, Dorothea; Eis-Hübinger, Anna M.; Kleim, Jörg-P.; Kaiser, Rolf; Bierhoff, Erhard; Schneeweis, Karl E. (1991): Diagnosis of acute and latent varicella-zoster virus infections using the polymerase chain reaction. In *Journal of medical virology* 35 (2), pp. 136–141.
3. Enders, G. (1984): Varicella-zoster virus infection in pregnancy. In *Progress in medical virology. Fortschritte der medizinischen Virusforschung. Progres en virologie medicale* 29, pp. 166–196.
4. Harper, David R.; Kangro, Hillar O.; Heath, Raymond B. (1988): Serological responses in varicella and zoster assayed by immunoblotting. In *Journal of medical virology* 25 (4), pp. 387–398.
5. Heath, R. B.; Kangro, H. O. (1990): Varizella zoster. In Arie J. Zuckerman (Ed.): *Principles and practice of clinical virology*. 2nd ed. Chichester: Wiley.
6. Hudnall, S. David; Stanberry, Lawrence R. (2006): Human Herpesvirus Infections. In Richard L. Guerrant, David H. Walker, Peter F. Weller (Eds.): *Tropical infectious diseases. Principles, pathogens & practice*. 2nd ed. Philadelphia: Churchill Livingstone, pp. 590–620.

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EIASON[®] VZV IgG

IVD

Enzyme immunoassay for the qualitative determination of IgG-class antibodies against Varicella-Zoster Virus (VZV) in human serum or plasma (citrate, heparin).


Instruction for use

For in-vitro use only


Product of

IASON GmbH
Feldkirchner Straße 4
A – 8054 Graz-Seiersberg
Tel.: +43 (0)316 28 43 00
Fax: +43 (0)316 28 43 00-113
Email: order@iason.eu
www.iason.eu

Used IFU symbols

Symbol	English	Symbol	English
	In vitro diagnostic medical device		Conjugate
	Order number		Cut Off Control
	Manufactured by		Negative Control
	Storage at		Positive Control
	European Conformity		Washing buffer
	Expiry date		Substrate solution
	Lot number		Stop solution
	Microplate		Sample diluent

Intended use

For in-vitro use only.

The EIASON® VZV IgG ELISA is intended for the qualitative determination of IgG class antibodies against Varicella-Zoster Virus (VZV) in human serum or plasma (citrate, heparin).

Introduction

Varicella-Zoster Virus (human herpes virus 3, HHV-3) belongs to the α -subfamily of herpesviridae. The virus particles measure about 145 nm in diameter. They consist of doublestranded DNA, are surrounded by an icosahedral protein capsid and an envelope which contains both host cells and viral components. The virus is usually transmitted in respiratory secretions, and a single serotype causes varicella (Chickenpox), a highly infectious childhood disease, and zoster (shingles), a neurodermic disease; both diseases are found worldwide. Varicella is the acute disease which follows primary contact with the virus, whereas zoster is the response of the partially immune host to a reactivation of the varicella virus present in the body in latent form.

Varicella is endemic, most commonly affected are children between 2 and 6 years of age. The course of disease is usually mild and complicated only in immunocompromised children. Rare fatal cases show multiple necrotic lesions in brain, lung (varicella pneumonia), kidneys (hemorrhagic nephritis), spleen, bone marrow, and occasionally in the intestinal tract. The lethality of varicella is below 0.1%. In the infrequent adult infections the disease is more severe, and complications are to be expected in about 5% of all cases. Zoster is of low incidence and appears with increasing frequency and severity with advancing age.

Usually the process remains localized, generalization is frequently encountered in a state of immunosuppression. Fatal cases are very rare and nearly always caused by an underlying disease.

Species	Disease	Symptoms	Transmission route
Varicella-Zoster Virus (VZV)	Varicella (primary infection)	Fever, itching exanthema and vesicular eruptions. Complication: pneumonia, inflammations of the brain or bacterial infections. VZV can cause malformations in infants and abortions if the mother has a primary infection in the first trimester of pregnancy.	Transmission by droplet; Vertically transmitted infection
	Herpes Zoster (Reactivation)	Fever, malaise, and gastrointestinal disturbance, painful cutaneous vesicular eruption, rash usual unilateral. Complication: Postherpetic neuralgia, infection of the eye, ear, and throat may occur if ophthalmic, trigeminal, or geniculate ganglia are involved; Ramsay-Hunt syndrome.	

The presence of pathogen or infection may be identified by:

- Microscopy
- PCR
- Serology: e.g. by ELISA

Assay principle

The qualitative immunoenzymatic determination of IgG-class antibodies against VZV is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a second washing step unbound conjugate is removed. The immune complex formed by the

bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA microwell plate reader.

Warnings and precautions

The EIASON® VZV IgG kit is for in vitro diagnostic use only and is not for internal use in humans or animals. This product must be used strictly in accordance with the instructions set out in the Package Insert. IASON will not be held responsible for any loss or damage (except as required by statute) caused, arising out of non-compliance with the instructions provided. CAUTION: this kit contains material of human and/or animal origin. Handle kit reagents as if capable of transmitting an infectious agent.

Source material from Human origin which is used in this kit was tested and found negative for HbsAG and HIV as well as for HCV antibodies. However, since there is no diagnostic procedure that excludes these agents with 100 percent certainty all components should be handled as potentially hazardous material.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. Disposal of kit reagents should be in accordance with local regulations.

Damaged test kit

In case of serious damage to the test kit or components, the company IASON must be notified in writing at least one week after receiving the kit. Severely damaged single components should not be used for the test run. They have to be stored until a final solution has been found.

Shelf life and storage of reagents

This kit is stable until the stated expiry date if stored as specified. Upon receipt, store all reagents at 2-8°C.

Specimen collection and preparation

Use human serum or plasma (citrate, heparin) samples with this assay. If the assay is performed within 5 days after sample collection, the specimen should be kept at 2-8°C; otherwise they should be aliquoted and stored deep-frozen (-70 to -20°C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing. Heat inactivation of samples is not recommended.

Sample dilution

Before assaying, all samples should be diluted 1+100 with **DIL**. Dispense 10 µL sample and 1 mL **DIL** into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

Materials provided

It is very important to bring all reagents, samples and controls (1-9) to room temperature (20-25°C) before starting the test run!

1. **MPL Varicella-Zoster Virus (VZV) Coated Microplate (IgG):** 12 breakapart 8-well snap-off strips coated with Varicella-Zoster Virus (VZV) antigen; in resealable

- aluminium foil. Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2-8 °C;
2. **DIL IgG Sample Diluent:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap.
 3. **STOP Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
 4. **WASH Washing buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap. Dilute **WASH** 1 + 19; e. g. 10 mL Washing Buffer + 190 mL distilled water. The diluted buffer is stable for 5 days at room temperature (20 - 25 °C). In case crystals appear in the concentrate, warm up the solution to 37 °C e.g. in a water bath. Mix well before dilution.
 5. **CONJ Varicella-Zoster Virus (VZV) anti-IgG Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgG in phosphate buffer (10 mM); coloured blue; ready to use; black cap.
 6. **SUB TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB); ready to use; yellow cap; < 5 % NMP. The reagent is ready to use and has to be stored at 2 - 8 °C, away from the light. The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.
 7. **PC Varicella-Zoster Virus (VZV) IgG Positive Control:** 1 bottle containing 2 mL control (human serum or plasma); coloured yellow; ready to use; red cap.
 8. **CUTOFF Varicella-Zoster Virus (VZV) IgG Cut-off Control:** 1 bottle containing 3 mL control (human serum or plasma); coloured yellow; ready to use; green cap.
 9. **NC Varicella-Zoster Virus (VZV) Negative Control:** 1 bottle containing 2 mL control (human serum or plasma); coloured yellow; ready to use; blue cap.

For potential hazardous substances please check the safety data sheet.

Materials supplied

- 1 cover foil
- 1 instruction for use (IFU)
- 1 plate layout

Materials required but not provided in the kit

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450/620nm
- Incubator 37°C
- Manual or automatic equipment for rinsing wells
- Pipettes to deliver volumes between 10 and 1000 µL
- Vortex tube mixer
- Deionised or (freshly) distilled water
- Disposable tubes
- Pipe stand
- Timer

Assay procedure

General remarks

1. All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
2. Once the test has been started, all steps should be completed without interruption.
3. Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination
4. Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
5. As a general rule the enzymatic reaction is linearly proportional to time and temperature.

Test preparation

Please read the test protocol carefully **before** performing the assay. Result reliability depends on strict adherence to the test protocol as described.

If performing the test on ELISA automatic systems we recommend to increase the washing steps from three to five and the volume of washing solution from 300 µL to 350 µL to avoid washing effects.

Prior to commencing the assay, the distribution and identification plan for all specimens and controls should be carefully established on the result sheet supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Adjust the incubator to 37° ± 1°C.

1. Dispense 100 µL [NC], [CUTOFF], [PC] and diluted samples into their respective wells. Leave well A1 for substrate blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour ± 5 min at 37 ± 1°C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of [WASH]. Avoid overflows from the reaction wells. The soak time between each wash cycle should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.
5. Dispense 100 µL [CONJ] into all wells except for the blank well (e.g. A1). Cover with foil.
6. **Incubate for 30 min at room temperature.** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL [SUB] into all wells.
9. **Incubate for exactly 15 min at room temperature in the dark.**
10. Dispense 100 µL [STOP] into all wells in the same order and at the same rate as for the [SUB]. Any blue colour developed during the incubation turns into yellow.
11. Measure the absorbance of the specimen at 450/620nm within 30 min after addition of the [STOP].

or automated on:

- IASON® PersonalLab
- IASON® Quardette
- IASON® Gladiator

Measurement

Adjust the ELISA Microwell Plate Reader **to zero** using the **substrate blank**.

If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for **NC**, **CUTOFF**, **PC** and patient sample in the plate layout.

Dual wavelength reading using 620 nm as reference wavelength is recommended.

Where applicable calculate the **mean absorbance values** of all duplicates.

Results

Run validation criteria

In order for an assay to be considered valid, the following criteria must be met:

- **Substrate blank:** Absorbance value < **0.100**
- **NC:** Absorbance value < **0.200** and **CUTOFF**
- **CUTOFF:** Absorbance value **0.150 – 1.300**
- **PC:** Absorbance value > **CUTOFF**

If these criteria are not met, the test is not valid and must be repeated.

Calculation of results

The cut-off is the mean absorbance value of the cut-off determinations.

Example: Absorbance value cut-off 0.44 + absorbance value cut-off 0.42 = 0.86 / 2 = 0.43
cut-off = 0.43

Results in units (IAU – IASON arbitrary units)

$$\frac{\text{Patient (mean) absorbance value} \times 10}{\text{cut-off}} = [\text{Units} = \text{IAU}]$$

Example:
$$\frac{1.591 \times 10}{0.43} = 37 \text{ IAU}$$

Interpretation of results

IAU	
Cut Off	10 -
Equivocal	9 - 11 Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks. If the result is equivocal again the sample is judged as negative .
Negative	< 9 The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.
Postiv	> 11 Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).

Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data.

In immunocompromised patients and newborns serological data only have restricted value.

Antibody isotypes and state of infection

Serology	Significance
IgM	Characteristic of the primary antibody response High IgM titer with low IgG titer: → suggests a current or very recent infection Rare: → persisting IgM
IgG	Characteristic of the secondary antibody response May persist for several years High IgG titer with low IgM titer: → may indicate a past infection
IgA	Produced in mucosal linings throughout the body (⇒ protective barrier) Usually produced early in the course of the infection

Quality control

The regular use of control samples at several analyte levels is advised to ensure day-to-day validity of results. Controls should be tested as unknowns. Quality Control charts should be maintained to follow the assay performance.

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid. In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or IASON GmbH directly.

Test characteristics

The results refer to the groups of samples investigated; these are not guaranteed specifications.

Precision

Intraassay Precision (n=24)			Interassay Precision (n=12)		
Probe	Mean [OD]	CV [%]	Probe	Mean [IAU]	CV [%]
1	0.437	4.55	1	24.63	1.85
2	1.014	3.71	2	25.75	6.81
3	1.071	3.17	3	5.26	5.35

Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 100.0% (95% confidence interval: 71.51% - 100.0%).

Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 100.0% (95% confidence interval: 98.5% - 100.0%).

Interferences

Interferences with hemolytic, lipemic or icteric sera are not observed up to a concentration of 10 mg/mL hemoglobin, 5 mg/mL triglycerides and 0.5 mg/mL bilirubin.

Cross reactivity

Investigation of a sample panel with antibody activities to potentially cross-reacting parameter did not reveal evidence of false-positive results due to cross-reactions.

Limitations of the procedure

Bacterial contamination or repeated freeze-thaw cycles of the specimen may affect the absorbance values.

Legal aspects

Reliability of results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact IASON.

Therapeutic consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under Reliability of Results. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point Therapeutic Consequences are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

Useful publications

1. Croen, Kenneth D.; Straus, Stephen E. (1991): Varicella-zoster virus latency. In *Annual review of microbiology* 45, pp. 265–282. DOI: 10.1146/annurev.mi.45.100191.001405.
2. Dlugosch, Dorothea; Eis-Hübinger, Anna M.; Kleim, Jörg-P.; Kaiser, Rolf; Bierhoff, Erhard; Schneweis, Karl E. (1991): Diagnosis of acute and latent varicella-zoster virus infections using the polymerase chain reaction. In *Journal of medical virology* 35 (2), pp. 136–141.
3. Enders, G. (1984): Varicella-zoster virus infection in pregnancy. In *Progress in medical virology. Fortschritte der medizinischen Virusforschung. Progres en virologie medicale* 29, pp. 166–196.
4. Harper, David R.; Kangro, Hillar O.; Heath, Raymond B. (1988): Serological responses in varicella and zoster assayed by immunoblotting. In *Journal of medical virology* 25 (4), pp. 387–398.
5. Heath, R. B.; Kangro, H. O. (1990): Varizella zoster. In Arie J. Zuckerman (Ed.): *Principles and practice of clinical virology*. 2nd ed. Chichester: Wiley.
6. Hudnall, S. David; Stanberry, Lawrence R. (2006): Human Herpesvirus Infections. In Richard L. Guerrant, David H. Walker, Peter F. Weller (Eds.): *Tropical infectious diseases. Principles, pathogens & practice*. 2nd ed. Philadelphia: Churchill Livingstone, pp. 590–620.

Pipetting scheme

Allow all reagents to reach room temperature before use

1. Predilution	1+100 Sample with DIL
2. Pipetting	Leave well A1 for blank NC, CUTOFF, PC, Sample 100µL
3. Incubation	1 hour at 37°C
4. Washing	wash 3 times: aspirate or decant add 300 µL* of diluted WASH repeat wash step 2 x and dry on absorbant material
5. Pipetting	CONJ except blank 100µL
6. Incubation	30 min at room temperature in the dark (20-25°C)
7. Washing	See point 4.
8. Pipetting	SUB 100µL
9. Incubation	15 min. at room temperature (20-25°C) in the dark
10. Pipetting	STOP 100µL
11. Reading	450nm (RF 620nm) Overrange Filter: 405nm, Factor 3 (depends on photometer), reading within 30 min. Calculation: $\frac{\text{Patient (mean) absorbance value} \times 10}{\text{Cut-off}} = [\text{Units} = \text{IAU}]$

*Please increase the volume of diluted **WASH** to 350 µL, if necessary

Expected values

	IAU
Cut Off	10
Equivocal	9 - 11
Negative	< 9
Postive	> 11