EIASON[®] Borrelia IgM

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Literatur

1. Asbrink, Eva; Hovmark, Anders (1988): Early and late cutaneous manifestations in Ixodes-borne borreliosis (erythema migrans borreliosis, Lyme borreliosis). In *Annals of the New York Academy of Sciences* 539, pp. 4–15.
2. Barbour, Alan G. (2006): Relapsing Fever and Other Borrelia Diseases. In Richard L. Guerrant, David H. Walker, Peter F. Weller (Eds.): *Tropical infectious diseases. Principles, pathogens & practice*. 2nd ed. Philadelphia: Churchill Livingstone, pp. 499–510.
3. Deutsche Borreliose-Gesellschaft e.V. (2011): Diagnostik und Therapie der Lyme-Borreliose. Leitlinien der Deutschen Borreliose-Gesellschaft.
4. Johnson, Russell C.; Schmid, George P.; Hyde, Fred W.; Steigerwalt, A. G.; Brenner, Don J. (1984): *Borrelia burgdorferi* sp. nov. Etiologic Agent of Lyme Disease. In *International journal of systematic bacteriology* 34 (4), pp. 496–497. DOI: 10.1099/00207713-34-4-496.
5. Karlsson, Mats; Hovind-Hougen, Kari; Svenungsson, Bo; Stiernstedt, Göran (1990): Cultivation and characterization of spirochetes from cerebrospinal fluid of patients with Lyme borreliosis. In *Journal of Clinical Microbiology* 28 (3), pp. 473–479.
6. Luther, Birgit; Moskophidis, Matthäus (1990): Antigenic cross-reactivity between *Borrelia burgdorferi*, *Borrelia recurrentis*, *Treponema pallidum*, and *Treponema phagedenis*. In *Zentralblatt für Bakteriologie : international journal of medical microbiology* 274 (2), pp. 214–226.
7. Magnarelli, Louis A.; Fikrig, Erol; Padula, Steven J.; Anderson, John F.; Flavell, Richard A. (1996): Use of recombinant antigens of *Borrelia burgdorferi* in serologic tests for diagnosis of Lyme borreliosis. In *Journal of Clinical Microbiology* 34 (2), pp. 237–240.
8. Preac-Mursic, V.; Wilske, Bettina; Schierz, G. (1986): European *Borrelia burgdorferi* isolated from humans and ticks culture conditions and antibiotic susceptibility. In *Zentralblatt für Bakteriologie, Mikrobiologie, und Hygiene. Series A, Medical microbiology, infectious diseases, virology, parasitology* 263 (1-2), pp. 112–118.
9. Wilske, Bettina; Preac-Mursic, V.; Fuchs, R.; Schierz, G. (1990): Diagnostik der Lyme Borreliose. Diagnose und Labor. In *Laboratoriumsblätter* 40, pp. 24–36.

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EIASON® Borrelia IgM

IVD

Enzyme immunoassay for the qualitative determination of IgM-class antibodies against *Borrelia burgdorferi* in human serum, plasma or CSF (cerebrospinal fluid)

Instruction for use

For in-vitro use only

Product of

IASON GmbH

Feldkirchner Straße 4

A – 8054 Graz-Seiersberg


Tel.: +43 (0)316 28 43 00

Fax: +43 (0)316 28 43 00-113


Email: order@iason.eu

www.iason.eu

REF E11-007-96

Used IFU symbols

Symbol	English	Symbol	English
	In vitro diagnostic device		Microplate
	Order number		Cut Off control
	Manufactured by		Negative control
	Storage at		Positive control
	EC Declaration of conformity		Wash solution
	Expiry date		Substrate
	Lot number		Stop solution
	Conjugate		Sample diluent

Intended Use

For in-vitro use only.

The EIASON® Borrelia IgM kit is intended for the qualitative determination of IgM class antibodies against *Borrelia burgdorferi* in human serum or plasma (citrate, heparin) and in cerebrospinal fluid (CSF) offering increased diagnostic specificity and sensitivity by employing immunodominant antigens. In combination with the recombinant IgG ELISA each clinical phase of Lyme disease - from Erythma chronicum migrans (ECM) over Bannwarth-Syndrom up to Lyme arthritis - may be detected even in borderland cases providing unerring results and diagnostic aid for the clinical staff.

In case of clinically suspected neuroborreliosis, the detection of intrathecally produced Borrelia-specific IgG and IgM antibodies is indicated. For this purpose serum-CSF pairs are analysed in parallel in the quantifiable EIASON® Borrelia burgdorferi IgG/ IgM ELISA. The result can then be confirmed by Blot.

First, serum and CSF are included in the standard Borrelia screening. If the CSF is positive, the values for total IgG and/or IgM and albumin must be requested for serum and CSF as the basis for the detection of the Borrelia specific antibody index.

By this method based on the Borrelia-specific antibody index, an intrathecally produced antibody fraction can be differentiated from a fraction that has passed the blood-CSF barrier passively to enter the subarachnoid space. Applying the limiting value QLim from the quotient diagram according to Reiber ensures that no false-negative antibody index values will be calculated in case of a polyspecific intrathecally immune reaction.

In cases of serum-CSF pairs showing a pathologically elevated antibody index in the screening test, confirmation with Blot is recommended.

Summary

Spirochetes are motile bacteria with a periplasmic axial filament. All pathogenic species belong to the family Treponemataceae, which includes the three genera: Treponema, Borrelia, and Leptospira. The Treponemae are extremely long, flexible, filamentous cells that are usually held in a characteristic spiral, or coiled-spring shape. Borreliae are the largest Treponemataceae with very coarse and irregular spirals. Borrelia burgdorferi, the causative agent of Lyme disease, is transmitted mainly by ticks but probably also by other blood-sucking insects. Habitats are the wooded, humid and temperate regions of North America, Europe, North Africa, Australia and Japan; foresters, farmers, and anybody entering infested forests may be affected.

The degree of contamination of ticks amounts to 3-60% dependent on seasonal and regional differences; up to 30% of the population may be infected (about 1500 cases annually in USA, several hundred in Europe).

Species	Diseases	Symptoms	Mechanism of infection
Borrelia burgdorferi	Lyme diseases	<p>Multisystemic illness</p> <p>Early manifestations: Erythema migrans; flu-like symptoms like fever, myalgias and headaches; carditis; meningitis</p> <p>Late manifestations: Chronic neurological, dermatological and / or cardiological symptoms; Acrodermatitis chronica atrophicans; Herxheimer reaction; Lyme-Arthritis</p>	Transmission by tick bites infected with Borrelia burgdorferi (Ixodes dammini in USA; Ixodes ricinus in Europe)

The presence of pathogen or infection may be identified by

- Microscopy
- PCR
- Serology: e.g. ELISA

Assay principle

The qualitative immunoenzymatic determination of specific antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a second washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA microwell plate reader.

Warnings and precautions

The EIASON® Borrelia IgM kit is for in vitro diagnostic use only and is not for internal use in humans or animals. This product must be used strictly in accordance with the instructions set out in the Package Insert. IASON will not be held responsible for any loss or damage (except as required by statute) caused, arising out of non-compliance with the instructions provided.

CAUTION: this kit contains material of human and/or animal origin. Handle kit reagents as if capable of transmitting an infectious agent.

Source material from Human origin which is used in this kit was tested and found negative for HbsAG and HIV as well as for HCV antibodies. However, since there is no diagnostic procedure that excludes these agents with 100 percent certainty all components should be handled as potentially hazardous material.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. Disposal of kit reagents should be in accordance with local regulations.

Damaged test kits

In case of serious damage to the test kit or components, the company IASON must be notified in writing at least one week after receiving the kit. Severely damaged single components should not be used for the test run. They must be stored until a final solution has been found. After that they should be disposed in accordance with official regulations.

Shelf life and storage of reagents

The reagents are stable up to the expiry date stated on the label when stored at 2 – 8 °C.

Specimen collection and preparation

Use human serum, plasma (citrate, heparin) and CSF (cerebrospinal fluid) samples with this assay. If the assay is performed within 5 days after sample collection, the samples should be kept at 2-8 °C; otherwise they should be aliquoted and stored deep-frozen (-70 – -20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing.

Heat inactivation of samples is not recommended.

Serum and CSF must be withdrawn from the patient at the same day. Parallel testing of serum and CSF in the same run is recommended. To increase diagnostic safety, patient sera and CSF should be tested in duplicate.

Sample dilution

Serum 1:101 Before assaying, all samples should be diluted 1+100 with **DIL**. Dispense 10 µL sample and 1 mL IgM **DIL** into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

CSF 1:2 50 µL of the green **DIL** are pipetted directly into the well. Then 50 µL CSF are added. Mix well! Please be aware of the correct order.

Materials

1. **MPL** Borrelia burgdorferi Coated Wells (IgM): 12 breakapart 8-well snap-off strips coated with Borrelia burgdorferi antigen; in resealable aluminium foil. The break-apart snap-off strips are coated with recombinant Borrelia burgdorferi antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2-8 °C.
2. **DIL** IgM Sample Diluent: 1 bottle containing 100 mL of ready to use buffer for sample dilution; pH 7.2 ± 0.2; coloured green; white cap.
3. **STOP** Stop Solution: 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L. Ready to use; red cap.
4. **WASH** Washing Buffer (20x conc.): 1 bottle containing 50 mL of a 20-fold concentrated buffer (pH 7.2 ± 0.2) for washing the wells; white cap. Dilute Washing Solution 1+19; e.g. 10 mL **WASH** + 190 ml fresh and germ free redistilled water. The diluted buffer is stable for 5 days at room temperature.
In case crystals appear in the concentrate, warm up the solution to 37 °C e.g. in a water bath. Mix well before dilution.
5. **CONJ** Borrelia burgdorferi anti-IgM Conjugate: 1 bottle containing 20 mL of peroxidase labelled antibody to human IgM; coloured red, Ready to use; black cap.
6. **SUB** TMB Substrate Solution: 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1%; ready to use; yellow cap; < 5% NMP.
The reagent is ready to use and has to be stored at 2-8 °C, away from the light. The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.

7. **PC** Borrelia burgdorferi IgM Positive Control: 1 vial containing 2 mL control (human serum or plasma); coloured yellow; ready to use; red cap.
8. **CUTOFF** Borrelia burgdorferi IgM Cut-off Control: 1 vial containing 3 mL control (human serum or plasma); coloured yellow; ready to use; green cap.
9. **NC** Borrelia burgdorferi IgM Negative Control: 1 vial containing 2 mL control (human serum or plasma); coloured yellow; ready to use; blue cap.

Materials and equipment needed

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450/620nm
- Incubator 37°C
- Manual or automatic equipment for rinsing wells
- Pipettes to deliver volumes between 10 and 1000 µL
- Vortex tube mixer
- Deionised or (freshly) distilled water
- Disposable tubes
- Timer

Assay procedure

General remarks

1. All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
2. Once the test has been started, all steps should be completed without interruption.
3. Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination
4. Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
5. As a general rule the enzymatic reaction is linearly proportional to time and temperature.

Test preparation

Please read the test protocol carefully **before** performing the assay. Result reliability depends on strict adherence to the test protocol as described.

If performing the test on ELISA automatic systems we recommend to increase the washing steps from three to five and the volume of washing solution from 300 µL to 350 µL to avoid washing effects.

Prior to commencing the assay, the distribution and identification plan for all specimens and controls should be carefully established on the result sheet supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Adjust the incubator to 37° ± 1°C.

1. Dispense 100 µL [NC] [PC] [CUTOFF] and diluted samples into their respective wells. Leave well A1 for substrate blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour ± 5 min at 37±1°C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of [WASH]. Avoid overflows from the reaction wells. The soak time between each wash cycle should be >5sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.
5. Dispense 100 µL [CONJ] into all wells except for the blank well (e.g. A1). Cover with foil
6. **Incubate for 30 min at room temperature. Do not expose to direct sunlight.**
7. Repeat step 4.
8. Dispense 100 µL [SUB] into all wells
9. **Incubate for exactly 15 min at room temperature in the dark.**
10. Dispense 100 µL STOP into all wells in the same order and at the same rate as for the SUB, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance of the specimen at 450/620nm within 30 min after addition of the [STOP].

or automated on:

- **IASON® PersonalLab**
- **IASON® Quardette**
- **IASON® Gladiator**

Measurement

Adjust the ELISA Microwell Plate Reader to zero using the substrate blank.

If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank, subtract its absorbance value substrate blank from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at 450 nm and record the absorbance values for each control and patient sample in the plate layout.

Dual wavelength reading using 620 nm as reference wavelength is recommended.

Where applicable calculate the mean absorbance values of all duplicates.

Results

Assay validation criteria

In order for an assay to be considered valid, the following criteria must be met:

- Substrate blank: Absorbance value < 0.100
- [NC]: Absorbance value < 0.200 and < [CUTOFF]
- [CUTOFF]: Absorbance value 0.150 – 1.300
- [PC]: Absorbance value > [CUTOFF]

If these criteria are not met, the test is not valid and must be repeated.

Calculation of results

1. The **CUTOFF** is the mean absorbance value of the **CUTOFF** determinations.

Example: Absorbance value **CUTOFF** 0.44 + absorbance value **CUTOFF** 0.42 = 0.86
/ 2 = 0.43

$$\text{CUTOFF} = 0.43$$

Results in IASON arbitrary units (IAU)

$\frac{\text{Patient (mean) absorbance value} \times 10}{\text{CUTOFF}} = [\text{IASON arbitrary units} = \text{IAU}]$

Example: $\frac{1.591 \times 10}{0.43} = 37 \text{ IAU}$

Interpretation of results

Cut-off	10 IAU	
Positive	> 11 IAU	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).
Equivocal	9 – 11 IAU	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks. If the result is equivocal again the sample is judged as negative .
Negative	< 9 IAU	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.
Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data. In immunocompromised patients and newborns serological data only have restricted value.		

Antibody isotypes and state of infection

Serology	Significance
IgM	Characteristic of the primary antibody response High IgM titer with low IgG titer: → suggests a current or very recent infection Rare: → persisting IgM
IgG	Characteristic of the secondary antibody response May persist for several years High IgG titer with low IgM titer: → may indicate a past infection

Evaluation of CSF results

Calculation of threshold values for sera and CSF for the qualitative evaluation:

IgM serum: threshold value = OD $\boxed{\text{CUTOFF}}$ = 10 IAU

IgM CSF: threshold value = 0.5 x OD $\boxed{\text{CUTOFF}}$ = 5 IAU

Interpretation or results

IgM Results	Interpretation and further actions
OD _{Serum} 1:101 < threshold value serum or OD _{Serum} / OD $\boxed{\text{CUTOFF}}$ < 10 IAU and OD _{CSF} 1:2 < threshold value CSF or OD _{CSF} / OD $\boxed{\text{CUTOFF}}$ < 5 IAU	No specific IgM antibodies in serum detected. No specific IgM antibodies in CSF detected. No indication of neuroborreliosis No determination of antibody index (AI) IgM
OD _{Serum} 1:101 ≥ threshold value serum or OD _{Serum} / OD $\boxed{\text{CUTOFF}}$ ≥ 10 IAU and OD _{CSF} 1:2 < threshold value CSF or OD _{CSF} / OD $\boxed{\text{CUTOFF}}$ < 5 IAU	Detection of specific IgM antibodies in serum No specific IgM antibodies in CSF detected. Indication of a fresh Borrelia infection No indication of neuroborreliosis No determination of antibody index (AI) IgM
OD _{Serum} 1:101 < threshold value serum or OD _{Serum} / OD $\boxed{\text{CUTOFF}}$ < 10 IAU and OD _{CSF} 1:2 ≥ threshold value CSF or OD _{CSF} / OD $\boxed{\text{CUTOFF}}$ ≥ 5 IAU	No specific IgM antibodies in serum detected. Detection of specific IgM antibodies in CSF Neuroborreliosis is possible Determination of antibody index (AI) IgM
OD _{Serum} 1:101 ≥ threshold value serum or OD _{Serum} / OD $\boxed{\text{CUTOFF}}$ ≥ 10 IAU and OD _{CSF} 1:2 ≥ threshold value CSF or OD _{CSF} / OD $\boxed{\text{CUTOFF}}$ ≥ 5 IAU	Detection of specific IgM antibodies in serum Detection of specific IgM antibodies in CSF Neuroborreliosis is possible Determination of antibody index (AI) IgM

The test can be terminated if serum and CSF give negative IgG and IgM results in screening. The Borrelia specific antibody index cannot be determined.

In case the CSF OD exceeds the threshold value, the procedure continues. In the following, IgG and IgM tests will be described in parallel. However, the antibody index is to be determined only for the immunoglobulin class in which a positive screening result for the CSF has been obtained.

Calculation of the specific IgG in IAU/mL

For the conversion of the ODs of CSF and serum into units, extinction values ≤ 2.500 are required.

If the OD value for the serum or CSF is larger than 2.500 or even outside the measurable range, the serum-CSF pair must be retested together in another run after appropriately higher dilution. The following dilutions are recommended:

OD	Serum dilution	CSF dilution (Preparation in the well)
Standard dilution	1:101	1:2 50 μ l DIL + 50 μ L CSF

Result of standard dilution (OD):	Recommendation for further dilutions:	CSF dilution (Preparation in the well)
OD ≤ 2.500	No further dilution	
2.500 < OD ≤ 3.000	1:500	1:5 80 μ L DIL + 20 μ L CSF
OD > 3.000	1:1000; 1:2000; 1:4000	Dilution in a polypropylen (PP) tube 1:10; 1:20; 1:40

The anti-Borrelia specific IgM has to be calculated from OD values of serum and CSF ≤ 2.500 .

If several values are in this range, the value with the OD closest to 1.000 is selected.

Note: As the cut-off in this formula the value of the **CUTOFF** (cut-off for serum) of the corresponding IgM determination is to be used for serum and for CSF. The threshold value for CSF from the qualitative evaluation is not valid here.

$$(1.1) \text{ spec. IgM serum [IAU/ml]} = \text{OD}_{\text{Serum}} \times 10 \times \text{dilution factor} / \text{cut-off} \times 101$$

$$(1.2) \text{ spec. IgM CSF [IAU/ml]} = \text{OD}_{\text{CSF}} \times 10 \times \text{dilution factor} / \text{cut-off} \times 101$$

Calculation of the specific IgM - quotient (Qspec IgM)

$$(2) \quad Q_{\text{spec}} (\text{IgM}) = \text{spec. IgM CSF [IAU/mL]} / \text{spec. IgM serum [IAU/mL]}$$

Calculation of the quotients QIgM, QAlb

Total IgM and albumin in both serum and CSF are determined by using a suitable method.

From the values obtained, the IgM quotient QIgM and the albumin quotient QAlb are calculated according to the following formulas:

$$(3) \quad \text{IgM-quotient } \mathbf{QIgM} = \text{IgM conc. CSF} / \text{IgM conc. serum}$$

$$(4) \quad \text{albumin quotient } \mathbf{QAlb} = \text{albumin conc. CSF} / \text{albumin conc. serum}$$

Calculation of the limiting QLim (IgM)

$$(5) \quad \mathbf{QLim} (\text{IgM}) = 0.67 \times \sqrt{\text{QAlb}^2 + (120 \times 10^{-6})} - 7.1 \times 10^{-3}$$

Calculation of the antibody index (AI)

The values obtained for QIgM and QLim (IgM) are compared and AI (IgM) is calculated using formula as follows:

(6) if $Q_{IgM} < Q_{Lim} (IgM)$: **AI (IgM)** = $Q_{spec} (IgM) / Q_{IgM}$

(7) If $Q_{IgM} > Q_{Lim} (IgM)$: **AI (IgM)** = $Q_{spec} (IgM) / Q_{Lim} (IgM)$

Interpretation of AI values

$0.6 < AI \leq 1.3$ (IgM) normal

$1.3 < AI < 1.5$ (IgM) borderline

$AI (IgM) \geq 1.5$ pathological

Antibody index values < 0.6 indicate an analytical fault (e.g. clinical chemistry). It is recommended to check whether serum and CSF are taken from the same day and if sample-mix-up can be excluded. A repetitive determination of the values of clinical chemistry and/ or of the Borrelia specific antibodies might be advisable.

Antibody index values $0.6 < AI < 1.3$ are in the normal range. An intrathecal synthesis of antibodies against Borrelia is unlikely. The detected Borrelia-antibodies in the CSF originate from the serum and have passively passed into the CSF space through the blood-liquor barrier.

Antibody index values between 1.3 and 1.5 are borderline cases. The serum-liquor pair should be retested. Or if clinical doubt exists, it is recommended to test another serum-CSF pair in the course of the infection.

Antibody index values ≥ 1.5 are to be classified as pathological. They indicate intrathecal synthesis of Borrelia-antibodies as found in patients suffering from neuroborreliosis.

Quality control

The regular use of control samples at several analyte levels is advised to ensure day-to-day validity of results. Controls should be tested as unknowns. Quality Control charts should be maintained to follow the assay performance.

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid. In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or IASON GmbH directly.

Test characteristics

Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte.

It is > 90 % for CSF samples and 98.35% (95% confidence interval: 95.84% - 99.55%) for serum samples.

Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is > 90 % for CSF samples and 92.45% (95% confidence interval: 81.79% - 97.91%) for serum samples.

Precision

Intra-Assay(n=24)			Inter-Assay (n=12)		
Sample	Mean (OD)	CV %	Sample	Mean (IAU)	CV %
1	0.451	4.58	1	26.33	11.26
2	1.172	7.76	2	39.3	8.31
3	1.840	10.22	3	6.38	12.82

Interferences

Interferences with hemolytic, lipemic or icteric sera are not observed up to a concentration of 10 mg/mL hemoglobin, 5 mg/mL triglycerides and 0.5 mg/mL bilirubin.

Note: The results refer to the groups of samples investigated; these are not guaranteed specifications.

Limitation of the procedure

Bacterial contamination or repeated freeze-thaw cycles of the specimen may affect the absorbance values.

By the use of recombinant antigens, cross reactions with antibodies against the following pathogens can be excluded, as best as possible:

- Treponema pallidum
- Leptospira
- Borrelia recurrentis

A lues-infection should be excluded, because antibodies against p41i could appear.

Legal aspects

Reliability of results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact IASON.

Therapeutic consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under Reliability of Results. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point Therapeutic Consequences are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

Useful publications

1. Asbrink, Eva; Hovmark, Anders (1988): Early and late cutaneous manifestations in Ixodes-borne borreliosis (erythema migrans borreliosis, Lyme borreliosis). In *Annals of the New York Academy of Sciences* 539, pp. 4–15.
2. Barbour, Alan G. (2006): *Relapsing Fever and Other Borrelia Diseases*. In Richard L. Guerrant, David H. Walker, Peter F. Weller (Eds.): *Tropical infectious diseases. Principles, pathogens & practice*. 2nd ed. Philadelphia: Churchill Livingstone, pp. 499–510.
3. Deutsche Borreliose-Gesellschaft e.V. (2011): *Diagnostik und Therapie der Lyme-Borreliose. Leitlinien der Deutschen Borreliose-Gesellschaft*.
4. Johnson, Russell C.; Schmid, George P.; Hyde, Fred W.; Steigerwalt, A. G.; Brenner, Don J. (1984): *Borrelia burgdorferi* sp. nov. Etiologic Agent of Lyme Disease. In *International journal of systematic bacteriology* 34 (4), pp. 496–497. DOI: 10.1099/00207713-34-4-496.
5. Karlsson, Mats; Hovind-Hougen, Kari; Svenungsson, Bo; Stiernstedt, Göran (1990): Cultivation and characterization of spirochetes from cerebrospinal fluid of patients with Lyme borreliosis. In *Journal of Clinical Microbiology* 28 (3), pp. 473–479.
6. Luther, Birgit; Moskophidis, Matthäus (1990): Antigenic cross-reactivity between *Borrelia burgdorferi*, *Borrelia recurrentis*, *Treponema pallidum*, and *Treponema phagedenis*. In *Zentralblatt für Bakteriologie : international journal of medical microbiology* 274 (2), pp. 214–226.
7. Magnarelli, Louis A.; Fikrig, Erol; Padula, Steven J.; Anderson, John F.; Flavell, Richard A. (1996): Use of recombinant antigens of *Borrelia burgdorferi* in serologic tests for diagnosis of Lyme borreliosis. In *Journal of Clinical Microbiology* 34 (2), pp. 237–240.
8. Preac-Mursic, V.; Wilske, Bettina; Schierz, G. (1986): European *Borrelia burgdorferi* isolated from humans and ticks culture conditions and antibiotic susceptibility. In *Zentralblatt für Bakteriologie, Mikrobiologie, und Hygiene. Series A, Medical microbiology, infectious diseases, virology, parasitology* 263 (1-2), pp. 112–118.
9. Wilske, Bettina; Preac-Mursic, V.; Fuchs, R.; Schierz, G. (1990): Diagnostik der Lyme Borreliose. Diagnose und Labor. In *Laboratoriumsblätter* 40, pp. 24–36.

Pipetting scheme

Allow all reagents to reach room temperature before use.

1. Predilution	1 part serum or plasma sample + 100 parts DIL	50µL DIL + 50 µL CSF
2. Pipetting	Leave well A1 for blank NC PC CUTOFF diluted sample 100µL	
3. Incubation	1 hour at 37°C	
4. Washing	wash 3 times: aspirate or decant add 300 µL * WASH repeat wash step 2 x and dry on absorbant material	
5. Pipetting	CONJ except blank	100µL
6. Incubation	30 min at room temperature (20-25°C)	
7. Washing	See point 4.	
8. Pipetting	SUB	100µL
9. Incubation	15 min. at room temperature (20-25°C) in the dark	
10. Pipetting	STOP	100µL
11. Reading	450nm (RF 620nm) Overrange Filter: 405nm, Factor 3 (depends on photometer), reading within 30 min. Calculation: Patient (mean) absorbance value x 10 / CUTOFF = [Units = IAU]	

*Please increase the volume of diluted WASH to 350 µL, if necessary

Expected values

Serum and plasma samples		Antibody Index IgM*	
	IAU		
Cut Off	10	Normal	0.6 – 1.3
Grey zone	9 - 11	Grey zone	1.3 – 1.5
Negative	< 9	Pathological	≥ 1.5
Postive	> 11		

*please refer to page 24 for details